

Ultra High Speed Imaging System Helps Unlock Understanding of Dynamics of Operation in Microdevices ...

Eberhard P. Hofer & Christian Rembe, Univ. of Ulm, Germany, and Joseph Honour, DRS Hadland Ltd., Tring, England

DRS Hadland, formerly Hadland Photonics, has supplied the Department of Measurement, Control and Microtechnology at the University of Ulm with an ultra-high speed digital imaging system to help them gain a greater understanding of the dynamics of moving parts in microdevices.

The measurement of position, velocity and acceleration of moving parts in microdevices is of great interest in the rapidly growing field of microsystem technology in order to solve the problems needed to allow evolution from laboratory use into industrial mass production. Such research is targeted towards gaining a better understanding of the dynamics of the moving parts in microdevices.

For the visual investigation of dynamics in microdevices high optical magnification is necessary because the dimensions of microsystem components are very small. Therefore, velocities of parts of such microdevices appear to be extremely fast and time constraints are very small.

The University of Ulm microtechnology research group has combined an intense pulsed light source, the ultra-high speed IMACON 468 camera from DRS Hadland and a standard optical microscope in a test rig to allow the analysis of nanosecond timescale dynamics of a series of different microdevices. To ensure optimal light input through the microscope a high pressure Xenon lamp was specially adapted to fit in its lamp


housing. The microscope was set up to image the magnified object onto the IMACON 468 ultra high speed camera. The prototype of this set-up for dynamic testing of microdevices was exhibited at the Hannover Industrial Fair in 1999.

The IMACON 468 received the Queens Award for Technology for its flexibility of operation and high performance in providing quantitative digital image analysis at high spatial resolution. Without mechanical moving parts the IMACON 468 ensured no vibration influences onto the tiny microdevices under study.

Using a proprietary design the IMACON 468 channels light through optical relays onto a special eight sided beamsplitter from where it passes onto eight intensified CCD units arranged in a circle around the beamsplitter. Production of eight discrete image frames through a single optical input not only provides very high spatial resolution but overcomes the drawback of image displacement optical anomalies caused by traditional high speed cameras use of several high speed camera mounted side-by-side.

Operating at up to 100 million frames per second using gated image intensifiers the IMACON 468 has provided the University of Ulm with a camera-microscope combination that is able to generate high contrast images (256 grey levels) with exposures as short as 10 ns and magnifications of up to 500 times.

Using this new technique Professor Hofer's research group has been able to analyse the motion of a 350 micron microturbine rotating at 200,000 revolutions per minute. In addition, they have looked at the heat induced stretch in a thermal micropump with a heating cycle of 2 microseconds, the oscillation of a microrelay with a switching time of 400 microseconds and to build up a detailed understanding of the mechanism of bubble-jet formation in thermal ink-jet printheads measuring 60 micron square. Using the technique to study microactuators, tiny devices that convert elec-



Launching a new product can be a treacherous journey.

Inspiration is not enough.

Success requires sound planning, experience, and great execution. If you are wise enough to recognize your organization's weaknesses, but not ready to hire full-time specialists, The Ascend Group may be your answer. We are a global network of professionals offering Sales, Marketing, and Strategic Management services to start-up and emerging technology companies.

Contact us today to find out about our range of available services.

Ascend
The Ascend Group
(503) 690-3555
Fax: (503) 690-1977
e-mail: Ascend@home.com

Taking your Business to the Top

trical energy into mechanical energy, the group has helped develop a new generation of diamond microheaters.

The IMACON 468 system has permitted previously impossible recording and analysis of the dynamics of operation of a series of microdevices by the Microtechnology group at Ulm which will lead to improvement of their design and, consequently, to the improvement of quality and lifetime of each microdevice. ■

Laboratory Manager Scanning Electron Microscopy Facility Smithsonian Institution

The National Museum of Natural History in Washington, DC is seeking an experienced electron microscopist to fill a vacancy for SEM laboratory operation and management. The SEM facility is designed to serve both the biological and geological research communities in the museum, and houses two recent model SEMs and one state-of-the-art environmental microscope. The principal responsibilities include training staff members and visiting scientists in proper use of equipment and theory of electron generation and detection, maintenance and troubleshooting all instrumentation, evaluation of new developments in SEM technology, and supervision of a support staff member. The successful applicant will also have the opportunity to gain experience in Focused Ion Beam (FIB) Microscopy.

This position will fill a federal government vacancy and is offered at the GS11 (\$42,724) or GS12 (\$51,204) levels, and carries promotion potential to grade GS 13. U.S. citizenship is required for this federal position. To obtain information concerning this vacancy call our automated Jobline at (202) 287-3102 (24 hours, 7 days a week), press 9 and request vacancy announcement **00MQ-2069**. Applications must be received by **May 16, 2000**. If questions arise after receiving and reading through the vacancy announcement please contact: Edward Vicenzi (Chair, SEM Lab Manager Search Committee) at *vicenzi@volcano.si.edu*. The Smithsonian Institution is an equal opportunity employer.



Smithsonian
National Museum of Natural History




**Problem: Find the right equipment,
... And find it now !!**

Solution: labgear.com

- ◆ Gather information about microscopy techniques
- ◆ Research and find analytical equipment
- ◆ Search by technique and then by manufacturer
- ◆ Video clips of laboratories and trade shows
- ◆ Calendar, new product features, tips from experts

20/20 Foresight Scientific Marketing
865/984-8683

NANO ANALYTICAL LABORATORY Freeze-Fracture Electron Microscopy of Pharmaceutical & Biological Samples



- *Including, but not limited to suspensions, gels or powders of
- *liposomes, proteo-liposomes, virosomes, niosomes, cochleate cylinders & lipids adapting non-bilayer structures,
- *applied to bacterial, cultured and/or tissue cells,
- *investigated at temperatures below/above lipid phase transition, and in media of any kind.

3951 Sacramento Street, San Francisco, CA 94118
TEL 415.666.3956 FAX 415.831.2813
E-mail *brigit@itsa.ucsf.edu*
<http://www.nanoanalytical.8m.com>