

# An Introduction to 3D Microscopy Techniques

Megan MacNeil and Duncan McMillan  
Carl Zeiss MicroImaging, Inc. Thornwood, NY.  
MMacNeil@zeiss.com

## From Confocal Microscopy to Deconvolution

In the earliest years of microscopy, it was well understood that visualization of the finest details was often limited by the 'intrusion' of light from out-of-focus planes. For many decades the only solution to this dilemma was to limit microscopic studies to very thin samples or to physically squash or slice thicker samples. More recently, biomedical research has driven the need for microscopes that can resolve very fine detail in 3 dimensions within intact, and often living, specimens. The use of fluorescence labeling further exacerbates the problem of out-of-focus light because signal is generated throughout the volume of the sample. The following article describes the most prevalent techniques for 3D imaging and important optical considerations for achieving highest image quality.

### Confocal

The distinctive optical design feature of a confocal laser scanning microscope (LSM) compared to a conventional microscope is the confocal aperture (or pinhole) arranged in a plane conjugate to the intermediate image plane and the object plane of the microscope. The detector only detects light that passes through the pinhole – *i.e.* light that emanates from the object plane in focus. Because light from above and below the object plane is blocked, the confocal microscope can be understood as an inherently depth-discriminating – or sectioning – optical system.

By varying the pinhole diameter, the thickness of the optical section can be adapted to practical requirements – for example, the pinhole diameter can be increased to collect a thicker section and hence more light. As an added advantage, the pinhole suppresses stray light, which improves image contrast. Application advantages of the confocal approach are the possibility to acquire multiple channels simultaneously (depending on the number of detectors), the use of laser for photo manipulation (*e.g.* FRET) and the immediate access to the confocal image. Disadvantages include the high light energy in the scanned pixel, the costs for lasers and the relative complexity of the systems

### 2-Photon

2-photon (2P) microscopy is founded on the fundamental principle of 2P excitation – that is, the absorption of two photons, at almost exactly the same time, with a fluorescent molecule. The combination of their energies causes the molecule to fluoresce as though it had been excited by a single photon of half the wavelength. At typical microscope light intensities the probability of this occurring is extremely low, but by temporally compressing the light into short, but very intense, pulses the photon density in the focal volume is high enough to generate significant fluorescence. However, outside of

the focal volume the photon density, and hence the generation of fluorescence, drops off dramatically. Therefore, while the 2P microscope is in many ways similar to a confocal microscope, the optical section formation itself is very different; in the confocal case fluorescence is generated throughout the vertical, or Z, axis and the pinhole eliminates the unwanted light from above and below the object plane. In the 2P case, fluorescence is *only* generated in the object plane, thereby eliminating the need for a pinhole.

This critical difference means that 2P imaging is a more efficient technique – only the fluorophores in the object plane are excited, thus reducing the net photo-bleaching and photo-toxicity caused during a 3D acquisition. The other key benefit of 2P imaging is deeper penetration (100's of microns) due to the reduced scattering of the longer wavelength (typically 700 – 1000 nm) light. 2P is the method of choice for imaging deep in opaque tissues, as necessary for intravital experiments.

### Structured illumination

In addition to the different laser-based optical sectioning techniques, there are also several widefield methods. One such technique is structured illumination or "grid projection." The principle of structured illumination centers around a grid of parallel lines projected into the image plane – thus, a high-resolution objective images the grid lines in the same focal plane as the sample. A CCD sensor captures this data as well as the out-of-focus data that stems from above or below the focal plane. To calculate the optical section, three images are captured, each with the grid shifted by one-third. Within this set of images, the out-of-focus data appears nearly identical. It is subtracted out of the image, whereas the in-focus data is an addition of the subtracted images.

Using this technique, there are several factors that determine the optical section thickness, including: magnification and numerical aperture of the objective, wavelength, and the line pairs per millimeter of the grid. The application limitations are the rela-




Image 1 (Widefield Image)




Image 2 (Deconvolved Image)




Image 3 (ApoTome section)





Image 4 (Confocal with Deconvolution)




*Caption: Images 1 to 4 show yeast cells with GFP transfected inner centromere proteins using conventional widefield (Image 1), Deconvolution in the Region of Interest (Image 2), Grating Imaging with the ApoTome (Image 3), and Confocal Imaging with Deconvolution in the Region of Interest (Image 4). Images 1 to 4 were acquired with the Plan-Apo 63x/1.4 NA Objective.*

# POWER THROUGH SIMPLICITY - AN IMAGING SOLUTION YOU COULD ONLY DREAM ABOUT



## ZEISS CENTRA® 100

An efficient, fully digital, and robust  
multi-purpose 100 kV TEM.  
Designed for rapid biomedical applications.  
CENTRA® 100 delivers outstanding image quality.



Enabling the Nano-Age World®

CENTRA® 100 – excellent image quality combined with ease of use and compact working area

Carl Zeiss SMT Inc.  
One Corporation Way  
Peabody, MA 01960  
USA

Tel. +1 978 / 826 7909  
Fax +1 978 / 532 5696  
info-usa@smt.zeiss.com  
www.smt.zeiss.com/nts



tive slow image acquisition and the need to acquire images sequentially (for fluorescence channels and grid positions). This limits the use for fast live cell imaging applications (>4fps), but allows an easy and convenient high resolution imaging of fixed or slow moving samples at very limited cost.

### Deconvolution

A computational method to optical sectioning is deconvolution. Unlike the methods described above, deconvolution is purely mathematical. The basis of the calculations is the point-spread function, or the diffraction pattern of a point of light in three dimensions. Point spread functions vary from system to system, as the properties of the objective, wavelength of light, alignment of the microscope all play a role. Blur characterized by the point spread function is mathematically modeled as a convolution, and the convolution is applied to every point of light within the image. Deconvolution reverses this phenomenon and deblurs or reconstructs the data.

There are many different published deconvolution algorithms, but they segment into two groups: deblurring or restorative. Deblurring methods are 2D methods, where each XY plane is deconvolved individually, using image information only from that plane or the closest neighboring planes. Deblurring calculations are much quicker calculations than their reconstructive counterparts. Restorative deconvolution, takes into account that blur is a 3D problem, and utilizes 3D data to subtract blur or reassign it to its original point of origin.

As an image processing approach, deconvolution does not immediately display a 'convolution free' image. But with increasing computer power this methods has again gained new application, especially for low light and fast imaging set ups. The strength of the restorative methods is to reassign the photons to the original source, gaining intensity and resolution.

### Important optical considerations for 3D applications

Microscope objectives are generally designed to ensure best optical quality just beneath the cover slip. To reach higher magnifications, an immersion media is applied, mostly oil with the refractive index of the glass. Imaging three dimensional and thicker specimens normally results in a gradual loss of image quality with increasing penetration depth into the sample. The signal to noise ratio decreases and the image is of lower contrast.

Other effects include the reduction of resolution as well as an increasing error in calibration precision in the Z-Axis. This results in a reduced Z-extension of the cell, if the image is displayed in a rendering program, so *e.g.* the nucleus is flattened in the Z-projection. All of these effects are caused by scattering light in the sample and on the surface borders between embedding media, cover slip and immersion media and is called spherical aberration. It can be easily observed by comparing the shape of the airy disc of a small fluorescence particle above and below the object. Spherical aberration is a problem if the disks are not symmetrical. In order to reduce this effect and to optimize image quality as well as image true spatial locations of the fluorophores different strategies can be applied:

Optimize the optical light path for minimal refractive index mismatch. The refractive index of the immersion media should resemble that of the embedding media as closely as possible. Most common are water- or glycerin immersion objectives.




Image 5 (Optic comparison oil-objective)




Image 6 (Optic comparison water-objective)

For a better adjustment, some of the objectives are equipped with a correction collar to precisely adjust the refractive index. Special live cell objectives combine high numerical apertures with Water or Glycerin Immersion (LCI Plan-Neofluar 63x/1.25 Imm. Corr.).

High aperture objectives are highly sensitive to differences in cover slip thickness. The use of closely tolerant cover slip thicknesses helps to avoid additional deteriorating effects in the images

Furthermore, imaging software can be used to correct for spherical aberration. The Deconvolution module in AxioVision offers the ability to correct for refractive Index mismatch, as well as the effects of increasing depth of the focal plane in the sample.

Please visit [www.zeiss.com/3dimaging](http://www.zeiss.com/3dimaging) to find out additional information on applications, techniques, and 3D imaging solutions from Carl Zeiss. Please visit [www.zeiss.com/objectives](http://www.zeiss.com/objectives) to for detailed information on microscope objectives from Carl Zeiss. ■

## SYMPOSIUM ANNOUNCEMENT

### Teaching Microscopy and Microanalysis

Microscopy & Microanalysis '08

August 3<sup>rd</sup> – 7<sup>th</sup>, Albuquerque, NM

- How do we learn microscopy and microanalysis methods?
- How do we effectively transfer our knowledge to those new to microscopy?
- How do we inspire students to enter and remain in the field?
- How do we keep up with advances in microscopy?

These are just a few of the questions that we hope to address in the symposium, "Teaching Microscopy and Microanalysis", to be held at M&M 2008 in Albuquerque, New Mexico. We will explore issues facing those interested in teaching and learning light, electron, and analytical microscopies. Best practices in classroom teaching, technician training, short courses, remote learning, community outreach and vendor training are among the areas that will be discussed. Contributed papers are welcome.

Submission and registration information are available in the M&M 2008 Call for Papers and on the Microscopy Society of America website, [www.microscopy.org](http://www.microscopy.org). Specific inquiries may be addressed to the organizers:


Charles Lyman, Lehigh University,  
[charles.lyman@lehigh.edu](mailto:charles.lyman@lehigh.edu)

Elaine Schumacher, McCrone Associates,  
[eschumacher@mccrone.com](mailto:eschumacher@mccrone.com)

# Detectors from Oxford Instruments - Performance beyond the expected

## Are you choosing the right EDS detector?

To ensure accurate and consistent results you can trust, you need the right detector for your application.



Email our specialists at [nanoanalysis@oxinst.com](mailto:nanoanalysis@oxinst.com) to find out which detector you should choose, or visit our website at [www.oxford-instruments.com/eds](http://www.oxford-instruments.com/eds)

The Business of Science™



**MICROSCOPY and MICROANALYSIS 2008**

**Mark Your Calendars!**

**ALBUQUERQUE**  
NEW MEXICO

August 3 - 7, 2008


### DID YOU RECEIVE TWO COPIES OF THIS ISSUE OF MICROSCOPY TODAY?

We have incorporated the USA attendance list from the Microscopy and Microanalysis-2007 meeting in Fort Lauderdale into the MT subscriber database. We have tried to eliminate duplicates, but we suspect that we were not 100% successful. If you received two copies of MT, please notify [microscopytoday@tampabay.rr.com](mailto:microscopytoday@tampabay.rr.com) to have the duplicate list entry removed. Thanks!

## MRS-5

We are ISO-9000 certified and ISO-17025 accredited  
**Microscopy Calibration Standard**  
Now you can calibrate from 1,000X to 1,000,000X!

This is our fourth generation, traceable, magnification reference standard for all types (SEM, FESEM, Optical, STM, AFM, etc.) of microscopy. The MRS-5 has multiple X and Y pitch patterns ranging from 80nm ( $\pm 1$ nm) to 2 $\mu$ m and 3 bar targets from 80nm to 3 $\mu$ m. There is also a STM test pattern.



Free web resource guide!



**GELLER  
MICROANALYTICAL  
LABORATORY, Inc.**

426e Boston St., Topsfield, Ma 01983  
[www.gellermicro.com](http://www.gellermicro.com)