

Complimentary XPS and AES Analysis of MoS₃ Solid Lubricant Coatings

S.S. Alnabulsi¹, J.R. Lince², D.F. Paul¹, J.F. Moulder¹, J.S. Hammond¹

¹ Physical Electronics, 18725 Lake Drive East, Chanhassen, MN 55317, USA

² Space Materials Laboratory, The Aerospace Corporation, El Segundo, CA 90245, USA

Molybdenum disulfide (MoS₂) nanoparticles are an ideal additive in solid coatings for lubricating mechanisms in vacuum environments, with widespread application in the spacecraft industry. The formation of these nanoparticles however can be complex, and the use of MoS₃ nanoparticles, which are produced using a simple wet chemical synthesis is being explored as an alternative.^[1] The use of MoS₃ as a tribological material has not been explored beyond its use as an oil additive.^[2] There is new interest in investigating its potential for use in solid lubricant coatings.

To aid in the evaluation of the tribological performance of a MoS₃-formulated coating compared to MoS₂ based coatings, x-ray photoelectron spectroscopy (XPS) and Auger electron spectroscopy (AES) are utilized as complimentary techniques for the surface characterization of the contact wear regions created on the coating surface.

The micro-XPS results provided quantitative chemical characterization that complement high spatial resolution AES imaging analysis of the molybdenum sulfide nanoparticles. A unique scanning micro-focused monochromatic x-ray source was used to provide x-ray excited secondary electron images (SXI) of the 50-100 μm wide wear track. The SXI images provide the ability to observe topography and regions with differences in photoelectron yield, making them very useful in locating features with different surface chemistry.

The MoS₃ nanoparticles were produced using simple wet chemical synthesis and are 20 nm to 30 nm in size. The nanoparticles are amorphous, and are contained in a dry coating comprised of a phenolic resin binder, similar to many commercial coatings used on spacecraft and in the aircraft industry.

The tribology of the resin-bonded MoS₃ nanoparticle coating is comparable to similarly prepared bonded coatings containing MoS₂. Tribometer testing showed the MoS₃-formulated coating performed similar to the MoS₂-based coatings, with similar coefficients of friction and endurances in dry nitrogen, with excellent solid lubricating properties.

The results of micro-area AES surface analyses within the contact wear region reveal changes in composition of the coating surface (Figure 1). The AES results are able to confirm the spatial distribution of MoS₃ particles in the binding material of the coating, and conversion into thin MoS₂ films in parts of the wear track.

The micro-area XPS surface analysis results of the wear region show a surface composition that is consistent with the production of a thin film of MoS₂ in the contact wear region (Figure 2). The results also show that after prolonged tribometer testing, the coating within the contact wear region does not fully convert to MoS₂.

References:

- [1] P. Afanasiev, "Synthetic approaches to the molybdenum sulfide materials," *Comptes Rendus Chimie*, 11(1–2) (2008) 159–182.
- [2] O.P. Parenago, V.N. Bakunin, G.N. Kuz'mina, A.Yu. Suslov, and L.M. Vedeneeva, "Molybdenum Sulfide Nanoparticles as New-Type Additives to Hydrocarbon Lubricants," *Doklady Chemistry*, 383(1–3) (2002) 86–88.




Figure 1. AES small area analysis within the wear track using a 3kV electron beam. The S map and spectral analysis highlight areas where MoS₂ has formed.




Figure 2. XPS Analysis of the track using a 20 μm scanning x-ray beam, shows an increase of MoS₂ formation in the contact region.