

Volume 23, Number 3 June 2017

# Microscopy AND Microanalysis



CAMBRIDGE  
UNIVERSITY PRESS

ISSN 1431-9276

# DIATOME

## DIAMOND KNIVES

# 40 YEARS

## OF DEVELOPMENT, MANUFACTURING, AND CUSTOMER SERVICE

What have we achieved in this period?

**ultra 45°** the first diamond knife with an absolutely score-free, hydrophilic cutting edge.

**semi** the first diamond knife for alternating sectioning ultrathin/semithin.

**cryo** the diamond knife for sectioning at low temperature.

**histo** the first diamond knife for semithin sections for light microscopy.

**ultra 35°** the diamond knife for optimized sectioning results in almost all applications.

**STATIC LINE II** the ionizer for eliminating electrostatic charging in ultramicrotomy.

**cryo-P** a cryo knife with a patented platform for section pick up.

**cryo immuno** the optimized cryo diamond knife for the Tokuyasu technique.

**ultra sonic** the oscillating diamond knife for room temperature sectioning.



**cryotrim** 45 and 25 optimizing trimming with diamond blades.

**ultra AFM & cryo AFM** the first diamond knives for AFM at room and low temperatures.

**cryo 25°** for sectioning frozen hydrated specimens.

What services can we offer you?

- Technical assistance in all fields of ultramicrotomy.
- Free sectioning tests for all types of samples.
- Make use of our many years of experience in perfecting our knives.
- Custom knives, tools, and boats.
- Special purchase programs.
- Workshops and training.



# DIATOME

for all your sectioning requirements

For more information, please call or write us today,  
or visit us online at: [www.emsdiasum.com](http://www.emsdiasum.com)

P.O. Box 550  
1560 Industry Rd.  
Hatfield, Pa 19440  
Tel: (215) 412-8390  
Fax: (215) 412-8450  
email: [sgkcck@aol.com](mailto:sgkcck@aol.com)  
[www.emsdiasum.com](http://www.emsdiasum.com)



# Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

THE OFFICIAL JOURNAL OF

MICROSCOPY SOCIETY OF AMERICA  
MICROANALYSIS SOCIETY  
MICROSCOPICAL SOCIETY OF CANADA /  
SOCIÉTÉ DE MICROSCOPIE DU CANADA  
MEXICAN MICROSCOPY SOCIETY  
BRAZILIAN SOCIETY FOR MICROSCOPY AND MICROANALYSIS  
VENEZUELAN SOCIETY OF ELECTRON MICROSCOPY  
EUROPEAN MICROBEAM ANALYSIS SOCIETY  
AUSTRALIAN MICROSCOPY AND MICROANALYSIS SOCIETY  
PORTUGUESE SOCIETY FOR MICROSCOPY

PUBLISHED IN AFFILIATION WITH

ROYAL MICROSCOPICAL SOCIETY  
GERMAN SOCIETY FOR ELECTRON MICROSCOPY  
BELGIAN SOCIETY FOR MICROSCOPY  
MICROSCOPY SOCIETY OF SOUTHERN AFRICA

## Editor-in-Chief

John Mansfield  
4304 Spring Lake Blvd.  
Ann Arbor, MI 48108-9657  
e-mail: jfmjfm@umich.edu

## Administrative Editor

John Shields  
University of Georgia  
Athens, GA 30602  
e-mail: jpshield@uga.edu

## Biological Sciences Applications Editors

W. Gray (Jay) Jerome  
Department of Pathology, Microbiology and  
Immunology  
U-2206 MCN  
Nashville, TN 37232-2561  
e-mail: jay.jerome@vanderbilt.edu

Heide Schatten  
Department of Veterinary Pathobiology  
University of Missouri-Columbia  
Columbia, MO 65211  
e-mail: SchattenH@missouri.edu

Rosemary White  
CSIRO Plant Industry  
Canberra, ACT 2601, Australia  
e-mail: Rosemary.white@csiro.au

Elizabeth Wright  
Department of Pediatrics  
School of Medicine  
Emory University  
Atlanta, GA 30322  
e-mail: erwright@emory.edu

## Cultural Heritage Applications Editor

Cathy Selvius DeRoo  
Detroit Institute of Arts  
5200 Woodward Avenue  
Detroit, Michigan 48202  
e-mail: cselviusderoo@gmail.com

## Materials Sciences Applications Editors

Vinayak Dravid  
Materials Science and Engineering  
Northwestern University  
Evanston, IL 60208-3105  
e-mail: v-dravid@northwestern.edu

Georg E. Fantner  
Interfaculty Institute for Bioengineering  
École Polytechnique Fédéral de Lausanne  
Lausanne, 1015 Switzerland.  
e-mail: georg.fantner@epfl.ch

David J. Larson  
CAMECA  
5500 Nobel Drive  
Madison, WI 53711  
e-mail: david.larson@ametec.com

Ross Marceau  
Institute for Frontier Materials  
Deakin University  
Geelong, VIC 3216, Australia  
e-mail: r.marceau@deakin.edu.au

Joseph Michael  
Sandia National Laboratories  
P.O. Box 5800  
Albuquerque, NM 87185  
e-mail: jrmicha@sandia.gov

Yoosuf N. Picard  
Materials Science & Engineering  
Carnegie Mellon University  
Pittsburgh, PA 15213  
e-mail: ypicard@cmu.edu

Masashi Watanabe  
Dept. of Mater. Sci. & Eng.  
Lehigh University  
Bethlehem, PA 18015  
e-mail: masashi.watanabe@lehigh.edu

## Special Issues and Reviews Editor

David J. Smith  
Department of Physics  
Arizona State University  
Tempe, AZ 85287-1504  
e-mail: david.smith@asu.edu

## Book Review Editor

Cynthia Goldsmith  
Centers for Disease Control  
Atlanta, GA 30333  
e-mail: csg1@cdc.gov

## M&M Program Guide Editor

Richard L. Martens  
1013 Bevell Building  
Box 870164  
Tuscaloosa, AL 35487-0164  
e-mail: rmartens@caf.ua.edu

## Proceedings Editor

Gail Celio  
University of Minnesota  
St. Paul, MN 55108  
e-mail: celio001@umn.edu

## Editorial Board

Ralph Albrecht	<i>University of Wisconsin, Madison, Wisconsin</i>
Ilke Arslan	<i>Pacific Northwest Laboratory, Richland, Washington</i>
Mary Grace Burke	<i>University of Manchester, Manchester, UK</i>
Barry Carter	<i>University of Connecticut, Storrs, Connecticut</i>
Wah Chiu	<i>Baylor College of Medicine, Houston, Texas</i>
Marc De Graef	<i>Carnegie Mellon University, Pittsburgh, Pennsylvania</i>
Niels de Jonge	<i>INM Institute for New Materials, Saarbrücken, Germany</i>
Elizabeth Dickey	<i>North Carolina State University, Raleigh</i>
Mark Ellisman	<i>University of California at San Diego, San Diego, California</i>
Pratibha Gai	<i>University of York, United Kingdom</i>
Marija Gajdardziska-Josifovska	<i>University of Wisconsin-Milwaukee, Milwaukee, Wisconsin</i>
Paul Kotula	<i>Sandia National Labs, Albuquerque, New Mexico</i>
William Landis	<i>University of Akron, Akron, Ohio</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Dale Newbury	<i>National Institute of Standards and Technology, Gaithersburg, Maryland</i>
Robert Price	<i>University of South Carolina, Columbia, South Carolina</i>
Jean-Paul Revel	<i>California Institute of Technology, Pasadena, California</i>
David Smith	<i>Arizona State University, Tempe, Arizona</i>
Nan Yao	<i>Princeton University, Princeton, New Jersey</i>
Nestor Zaluzec	<i>Argonne National Laboratory, Argonne, Illinois</i>

## Editorial Board Representatives from Affiliated Societies

Masashi Watanabe	<i>Lehigh University (MAS)</i>
Gautam Kumar Dey	<i>Bhabha Atomic Research Centre (EMSI)</i>
Gema Gonzalez	<i>Venezuelan Institute for Scientific Investigation (Venezuela)</i>
Michael Robertson	<i>Acadia University, Wolfville, Nova Scotia (Canada)</i>
Brendan Griffin	<i>University of Western Australia (AMMS)</i>
Guillermo Solorzano	<i>Pontificia Universidade Catolica, Rio de Janeiro (Brazil)</i>
Mike Matthews	<i>Atomic Weapons Establishment, Reading, Great Britain (EMAS)</i>
Miguel Yacamán	<i>Mexico Institute for Nuclear Research (Mexico)</i>
Henrique Almeida	<i>Universidade do Porto (Portugal)</i>

## Founding Editor

Jean-Paul Revel	<i>California Institute of Technology, Pasadena, California</i>
-----------------	---

## Previous Editors-in-Chief

Dale Johnson	<i>University of South Florida, Tampa, Florida</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Robert L. Price	<i>University of South Carolina, Columbia, South Carolina</i>

This journal is part of the **Cambridge Core** service. Access to online tables of contents and article abstracts is available to all researchers at no cost. Access to full-text articles online is provided to those with online subscription. Online subscriptions must be activated. Once your subscription is activated, free access to past, present, and forthcoming articles is available at:


**Microscopy and Microanalysis website: [cambridge.org/MAM](http://cambridge.org/MAM).**

Instructions for authors submitting manuscripts may be found at [cambridge.org/MAM](http://cambridge.org/MAM). Select "Further Information" then select "Instructions for Contributors." An abbreviated version of these instructions will be published in the first issue (February) of each volume.



# The NEXT Generation of 120 kV TEM

Innovation, Performance, and Reliability—The All New Hitachi *HT7800* Series



[Samples]

Top: Fin Field Effect Transistor  
Middle: Mouse Kidney  
Bottom: Single Crystal Silicon

Hitachi's advanced electron optics and dual-mode lens technology redefine the boundaries of 120 kV TEM performance. With  $\leq 0.19$  nm on-axis resolution, the highest contrast in its class, and diverse analytical capabilities for Material Science and Life Science fields, the *HT7800* Series provides a complete solution for all your TEM imaging needs.

For more details, contact us at: [microscopy@hitachi-hta.com](mailto:microscopy@hitachi-hta.com)

## Think Outside the Lab



*Microscopy and Microanalysis* publishes original research papers dealing with a broad range of topics in microscopy and microanalysis. These include articles describing new techniques or instrumentation and their applications, as well as papers in which established methods of microscopy or microanalysis are applied to important problems in the fields of biology or materials science. Microscopy and microanalysis are defined here in a broad sense, and include all current and developing approaches to the imaging and analysis of microstructure. The criteria for acceptance of manuscripts are the originality and significance of the research, the quality of the microscopy or microanalysis involved, and the interest for our readership.

Four types of communications are published in the Journal. **Regular Articles** are of substantial length and describe the findings of an original research project that satisfies the aims and scope of the Journal, described above. **Review Articles** summarize the current status of an important area within the aims and scope of the Journal. **Letters to the Editor** usually contain comments on recent articles that have appeared in the Journal. **Book Reviews** are also published, but these are solicited only through the Book Review Editor.

### Instructions for Contributors

Instructions for authors contributing manuscripts may be found at <http://mc.manuscriptcentral.com/mam> under "Resources: Instructions and Forms." Authors may also visit [cambridge.org/mam](http://cambridge.org/mam), select "Information," and then select "Instructions for Contributors." An abbreviated version of these instructions will be published in the first issue (February) of each volume.

### Copyright Information

Submission of a manuscript implies: that the work described has not been published before (except in the form of an abstract or as part of a published lecture, review, or thesis); that it is not under consideration for publication elsewhere; that its publication has been approved by all coauthors, if any, as well as by the responsible authorities at the institute where the work has been carried out; that, if and when the manuscript is accepted for publication, the authors agree to automatic transfer of the copyright to the Microscopy Society of America; that the manuscript will not be published elsewhere in any language without the consent of the copyright holders; and that written permission of the copyright holder is obtained by the authors for material used from other copyrighted sources.

All articles published in this journal are protected by copyright, which covers the exclusive rights to reproduce and distribute the article (e.g., as offprints), as well as all translation rights. No material published in this journal may be reproduced photographically or stored on microfilm, in electronic data bases, video disks, etc., without first obtaining written permission from the publisher.

The use of general descriptive names, trade names, trademarks, etc., in this publication, even if not specifically identified, does not imply that these names lack protection by the relevant laws and regulation.

Authorization to photocopy items for internal or personal use, or the internal or personal use of specific clients, is granted by Cambridge University Press, provided that the appropriate fee is paid directly to Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, USA (Tel: (508) 750-8400), stating the ISSN (1431-9276), the volume, and the first and last page numbers of each article copied. The copyright owner's consent does not include copying for general distribution, promotion, new works, or resale. In these cases, specific written permission must first be obtained from the publisher.

### Disclaimer

The Microscopy Society of America, the other societies stated, and Cambridge University Press cannot be held responsible for errors or for any consequences arising from the use of the information contained in this journal. The appearance of scientific reports and/or workshops, or any other material in *Microscopy and Microanalysis* does not constitute an endorsement or approval by The Microscopy Society of America of the findings, data, conclusions, recommendations, procedures, results, or any other aspect of the content of such articles. The appearance of advertising in *Microscopy and Microanalysis* does not constitute an

endorsement or approval by The Microscopy Society of America of the quality or value of the products advertised or any of the claims, data, conclusions, recommendations, procedures, results, or any other information included in the advertisements.

While the advice and information in this journal is believed to be true and accurate at the date of its going to press, neither the authors, the editors, nor the publisher can accept any legal responsibility for any errors or omissions that may be made.

### Subscription Information

*Microscopy and Microanalysis* is published bimonthly in February, April, June, August, October, and December by Cambridge University Press. Three supplements (*Meeting Guide*, *Program Guide*, and *Proceedings*) are published in June and August.

**Society Rates:** Members of the Microscopy Society of America should contact the MSA Business Office for all subscription inquiries: Microscopy Society of America, 11130 Sunrise Valley Dr, Suite 350, Reston, VA 20191, Tel.: (703) 234-4115, Email: [associationmanagement@microscopy.org](mailto:associationmanagement@microscopy.org), URL: [www.microscopy.org](http://www.microscopy.org). Members of other affiliated societies should contact their respective society business offices for all subscription inquiries.

**Subscription Rates:** Institutions print and electronic: US \$2251.00 in the USA, Canada, and Mexico; UK £1354.00 + VAT elsewhere. Institutions online only: US \$1460.00 in the USA, Canada, and Mexico; UK £883.00 + VAT elsewhere. Individuals print plus online: US \$658.00 in the USA, Canada, and Mexico; UK £400.00 + VAT elsewhere. Prices include postage and insurance.

**USA, Canada, and Mexico:** Subscribers in the USA, Canada, and Mexico should send their orders, with payment in US dollars or the equivalent value in Canadian dollars, to: Cambridge University Press, Customer Services Department (Journals), 1 Liberty Plaza, New York, NY 10006, USA. Tel: (845) 353-7500. Fax: (845) 353-4141. Orders may be phoned direct (toll free): (800) 872-7423. E-mail: [journals\\_subscriptions@cup.org](mailto:journals_subscriptions@cup.org).

**Outside North America:** Subscribers elsewhere should send their orders, with payment in sterling, to: Customer Services Department (Journals), Cambridge University Press, University Printing House, Shaftesbury Road, Cambridge, CB2 8BS, UK. Tel: +44 (0)1223 326070. Fax: 44 (0)1223 325150. E-mail: [journals@cambridge.org](mailto:journals@cambridge.org)

**Change of address:** Allow six weeks for all changes to become effective. All communications should include both old and new addresses (with postal codes) and should be accompanied by a mailing label from a recent issue. Society members should contact their respective society business offices to inform them of address changes.

### Editorial Office

John Mansfield, Editor in Chief, 4304 Spring Lake Blvd., Ann Arbor, MI 48108-9657, USA; Tel: (734) 936-3352; Fax: (734) 763-2282; E-mail: [jfmjfm@umich.edu](mailto:jfmjfm@umich.edu).

### Office of Publication


Cambridge University Press, 1 Liberty Plaza, New York, NY 10006, USA; Tel: (212) 337-5000; Fax: (212) 337-5959.

### Advertising Sales & Production

Kelly Miller, M.J. Mrvica Associates, Inc., 2 West Taunton Avenue, Berlin, NJ 08009, USA; Tel: (856) 768-9360; Fax: (856) 753-0064.

© 2017 by Microscopy Society of America. Printed in the United States on acid-free paper. Periodicals postage paid at New York, NY, and additional mailing offices. Return postage guaranteed. Postmaster: Send address changes in the U.S.A. and Canada to *Microscopy and Microanalysis*, Subscription Department, Cambridge University Press, 1 Liberty Plaza, New York, NY 10006.

# “For the best TKD/EBSD data, I always use **Evactron**<sup>®</sup> cleaning.”



88% increase in hit rate on 20nm Au with Bruker eFlash<sup>FS</sup> and OPTIMUS<sup>™</sup> TKD

Plasma cleaning improves TKD pattern quality which equals a better indexing rate, increased acquisition speed and improved effective spatial resolution. Evactron EP dual-action, *turbo plasma cleaning*<sup>™</sup> has:

- ◆ plasma etch plus UV active desorption
- ◆ quick pump down to high vacuum
- ◆ easy to use sample and chamber cleaning
- ◆ a compact and efficient remote plasma source

**Need clean sample surfaces?**

**Let us find the fastest way to pristine for you!**



**WWW.EVACTRON.COM**  
**1-650-369-0133**



## Precision, Speed, Stability

### NANO-POSITIONING FOR MICROSCOPY & IMAGING




Manual Stage  
25x25mm Travel




High-Speed Nanofocus  
<1nm Resolution



XYZ Flexure Stage  
<1nm Resolution




UHV Piezo Stage  
0.1nm Resolution, XYZ



XY Piezo Motor  
Stage, 100mm


**The Broadest and Deepest Portfolio of Precision Motion Technologies**



Vacuum Stages



Piezo Mechanics




Air Bearings



Fiber Alignment




Hexapods



Beam Steering

Learn more ► [www.pi-usa.us](http://www.pi-usa.us)  
508.832.3456  
PI (Physik Instrumente)





FEI Avizo® 3D visualization of two large adjacent crystalline dendrites of a bulk-metallic-glass matrix composite ( $Zr_{58.5}Ti_{14.3}Nb_{5.2}Cu_{6.1}Ni_{4.9}Be_{1.0}$ ). Data was obtained by large volume serial sectioning tomography using the Helios PFIB DualBeam. The sectioned block is about  $90 \times 80 \times 70 \mu m^3$ . Sample from The University of Tennessee, USA. Images courtesy of The University of Manchester.

## Large 3D volumes with unprecedented surface resolution

Until recently, available technologies have limited the volumes and depths of materials that can be analyzed at high resolution, ultimately restricting the insight into structural, crystallographic, and chemical properties. This is no longer the case. The Helios™ PFIB DualBeam™ offers unrivaled access to regions of interest deep below the surface—combining serial section tomography with statistically relevant data analysis.

Discover more at [FEI.com/Helios-PFIB](https://FEI.com/Helios-PFIB)



# Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

Volume 23, Number 3

June 2017

## MATERIALS SCIENCE APPLICATIONS

The Relationship Between Atomic Structure and Strain Distribution of Misfit Dislocation Cores at Cubic Heteroepitaxial Interfaces 449

*Cai Wen*

Influence of Noise-Generating Factors on Cross-Correlation Electron Backscatter Diffraction (EBSD) Measurement of Geometrically Necessary Dislocations (GNDs) 460

*Landon T. Hansen, Brian E. Jackson, David T. Fullwood, Stuart I. Wright, Marc De Graef, Eric R. Homer, and Robert H. Wagoner*

Method for Estimating the Charge Density Distribution on a Dielectric Surface 472

*Takuya Nakashima, Hiroyuki Suhara, Hidekazu Murata, and Hiroshi Shimoyama*

Mitigating Curtaining Artifacts During Ga FIB TEM Lamella Preparation of a 14 nm FinFET Device 484

*Andrey Denisjuk, Tomáš Hrnčíř, Jozef Vincenc Oboňa, Sharang, Martin Petrevec, and Jan Michalička*

Overcoming Peak Overlaps in Titanium- and Vanadium-Bearing Materials with Multiple Linear Least Squares Fitting 491

*Michael Mengason and Nicholas Ritchie*

Stability of a Bifunctional Cu-Based Core@Zeolite Shell Catalyst for Dimethyl Ether Synthesis Under Redox Conditions Studied by Environmental Transmission Electron Microscopy and *In Situ* X-Ray Ptychography 501

*Sina Baier, Christian D. Damsgaard, Michael Klumpp, Juliane Reinhardt, Thomas Sheppard, Zoltan Balogh, Takeshi Kasama, Federico Benzi, Jakob B. Wagner, Wilhelm Schwieger, Christian G. Schroer, and Jan-Dierk Grunwaldt*

Measurement Error in Atomic-Scale Scanning Transmission Electron Microscopy Energy-Dispersive X-Ray Spectroscopy (STEM-EDS) Mapping of a Model Oxide Interface 513

*Steven R. Spurgeon, Yingge Du, and Scott A. Chambers*

Local Structural Damage Evaluation of a C/C-SiC Ceramic Matrix Composite 518

*Benedicta D. Arhatari, Matthew Zonneveldt, John Thornton, and Brian Abbey*

Advanced Identification and Quantification of In-Bearing Minerals by Scanning Electron Microscope-Based Image Analysis 527

*Kai Bachmann, Max Frenzel, Joachim Krause, and Jens Gutzmer*

## BIOLOGICAL SCIENCE APPLICATIONS

The Application of Contrast Media for *In Vivo* Feature Enhancement in X-Ray Computed Tomography of Soil-Grown Plant Roots 538

*Samuel D. Keyes, Neil J. Gostling, Jessica H. Cheung, Tiina Roose, Ian Sinclair, and Alan Marchant*

New Methods of Morphometric Analyses on Scyphozoan Jellyfish Statoliths Including the first Direct Evidence for Statolith Growth Using Calcein as a Fluorescent Marker 553

*Ilka Sötje, Tamar Dishon, Frank Hoffmann, and Sabine Holst*

Segmentation and Quantitative Analysis of Apoptosis of Chinese Hamster Ovary Cells from Fluorescence Microscopy Images 569

*Yuncheng Du, Hector M. Budman, and Thomas A. Duever*

**Microscopy and Microanalysis website:** <http://www.journals.cambridge.org/MAM>  
Indexed in Chemical Abstracts, Current Contents, BIOSIS, and MEDLINE (PubMed)



**On the Cover:** Rouleaux formation, the aggregation of erythrocytes to resemble a stack of coins, is present when certain synthetic progestins used in combined oral contraceptives (COC) are added to whole blood. The cover image was visualised with a Zeiss Crossbeam 540 High Resolution SEM at 1 kV, and along with spontaneous fibrin formation also visualised with the same SEM may shed some light on the increased risk of venous thrombosis (VT) associated with COC use. For further details, see Swanepoel et al., pp. 607–617.

Local Response of Sialoliths to Lithotripsy: Cues on Fragmentation Outcome 584  
*Pedro Nolasco, Ana J. dos Anjos, José Dias, Paulo V. Coelho, Carla Coelho, Manuel Evaristo, Albano Cavaleiro, António Maurício, Manuel F. C. Pereira, Virgínia Infante, António P. Alves de Matos, Raúl C. Martins, and Patrícia A. Carvalho*

The Effect of Endogenous and Synthetic Estrogens on Whole Blood Clot Formation and Erythrocyte Structure 599  
*Albe C. Swanepoel, Odette Emmerson, and Ethersia Pretorius*

Effect of Progesterone and Synthetic Progestins on Whole Blood Clot Formation and Erythrocyte Structure 607  
*Albe C. Swanepoel, Odette Emmerson, and Ethersia Pretorius*

## INSTRUMENTATION AND SOFTWARE

Efficient Simulation of Secondary Fluorescence Via NIST DTSA-II Monte Carlo 618  
*Nicholas W. M. Ritchie*

PENEPMA: A Monte Carlo Program for the Simulation of X-Ray Emission in Electron Probe Microanalysis 634  
*Xavier Llovet and Francesc Salvat*

Pattern Center and Distortion Determined from Faint, Diffuse Electron Diffraction Rings from Amorphous Materials 647  
*János L. Lábár and Partha P. Das*

Measuring the Autocorrelation Function of Nanoscale Three-Dimensional Density Distribution in Individual Cells Using Scanning Transmission Electron Microscopy, Atomic Force Microscopy, and a New Deconvolution Algorithm 661  
*Yue Li, Di Zhang, Ilker Capoglu, Karl A. Hujsak, Dhwanil Damania, Lusik Cherkezyan, Eric Roth, Reiner Bleher, Jinsong S. Wu, Hariharan Subramanian, Vinayak P. Dravid, and Vadim Backman*

## MICROGRAPHIA


Effect of Sublethal Nickel Chloride Exposure on Crayfish, *Astacus leptodactylus* Ovary: An Ultrastructural, Autometallographic, and Electrophoretic Analyses 668  
*Otilia Zarnescu, Ana-Maria Petrescu, Alexandra Gaspar, and Oana Craciunescu*

Scanning Electron Microscopy Coupled with Energy Dispersive Spectrometric Analysis Reveals for the First Time Weddellite and Sylvite Crystals on the Surface of Involucral Bracts and Petals of two *Xeranthemum* L. (Compositae) Species 679  
*Milan Gavrilović, Suzana Erić, Petar D. Marin, Núria Garcia-Jacas, Alfonso Susanna, and Pedja Janačković*

Microstructure and Quantitative Micromechanical Analysis of Wood Cell-Emulsion Polymer Isocyanate and Urea-Formaldehyde Interphases 687  
*Lizhe Qin, Lanying Lin, Feng Fu, and Mizi Fan*



Someone has to be first.  
XFlash<sup>®</sup> FlatQUAD, multiple detector systems & VZ



Analyze textured  
samples with ease.




- The unique XFlash<sup>®</sup> FlatQUAD offers a large solid angle and an amazing take-off angle. This minimizes shadowing so you can see every nook and cranny of your sample in record time.
- Bruker's multiple detector systems provide large solid angles and great take-off angles too. You can start with one detector building up to four, depending on your needs.
- The Variable Z (VZ) adapter allows you to optimize take-off angles in-situ, which significantly improves the analysis of topographically challenging samples.

**Someone has to be first.**



[www.bruker.com/quantax-flatquad](http://www.bruker.com/quantax-flatquad)

# The Swiss Army Knife of Microscopy



Always have the right  
tool on hand!





Visit **TESCAN USA** at **M&M 2017**

[www.tescan-usa.com](http://www.tescan-usa.com) | 724-772-7433



# SIA

**1 TO 50 MEGAPIXELS** live and slow scan  
**MAGNIFICATION FACTOR OF 1** on bottom mounted cameras  
**DIFFRACTION BEAM STOP** on side mounted cameras



*Affordable* TEM camera systems for research, education, healthcare, and industry since 2001

**Scientific Instruments and Applications**  
2773 Heath Lane • Duluth, GA 30096  
(770) 232 7785 • [www.sia-cam.com](http://www.sia-cam.com)

## SEM & TEM Sample Preparation

### Microscopy



#### DESK V TSC Sample Coater

- Fully automated
- Minimize sample charge
- Reduce e-beam damage
- Increase secondary electron emissions

#### BENCH TOP TURBO Evaporation Tool for Tissue and Life Science

- Carbon coating
- Sputtering
- AC glow discharge
- Low-voltage thermal evaporation

### Support Performance Value

#### DV-502

Large Sample, High Throughput  
Sputter & Evaporation

- Rotary shadowing for metals and carbon
- High vacuum gate valve for rapid cycling
- Recipe development

**DENTON VACUUM**  
Enabling Innovation

All Denton Vacuum Microscopy tools are backed by a 12 month warranty, over 50 years of process knowledge, an in-house process engineering group, worldwide support, and Global factory service center.

To learn more about how to better prepare your samples, contact your electron microscope manufacturer or call us at +1-800-666-6004

[www.dentonvacuum.com](http://www.dentonvacuum.com)



# Deepen your knowledge.

ZEISS Crossbeam 550

// INNOVATION  
MADE BY ZEISS



## Your FIB-SEM for high throughput 3D analysis and nanofabrication

With ZEISS Crossbeam 550 you experience an outstanding combination of imaging performance and analytical power. Process any sample on a nanoscopic scale with the next-generation FIB. **Visit the ZEISS booth #1018 at M&M and see Crossbeam 550 in action.**

[www.zeiss.com/crossbeam](http://www.zeiss.com/crossbeam)



## IN SITU MADE SIMPLE


---

*In situ* microscopy just got a lot easier. Meticulously crafted by our world-class engineers, the next-generation Protochips holders deliver heating and electrical biasing in liquid or gas to your electron microscope with the click of a button. Self aligning parts and intuitive software controls mean you can focus on what matters, your experiment. Find out more at:

[www.protochips.com/products](http://www.protochips.com/products)

 **Protochips**  
Quantifiably Better™





## EDAX EDS Systems with New SDD Options Advanced Analysis with Superior Results

- Choice of optimized SDDs to suit your materials analysis needs
- Best light element sensitivity with silicon nitride ( $\text{Si}_3\text{N}_4$ ) window
- Vacuum encapsulated module
- Highest throughput SDD available, with unparalleled resolution
- Safe for plasma cleaning