










Powder Diffraction

An international
journal of materials
characterization


Volume 13 Number 4 December 1998





- optimum goniometer and sample accessibility and visibility
- stepper motors with optical encoders for highest precision
- dovetail-tracks for continuously adjustable measurement circle
- maximum modularity
- exchange of optical components in seconds
- exchange of detectors in seconds
- ceramic X-ray tube with reproducible focus position - 100% compatible with glass tubes
- Dynamic Scintillation Detector with large dynamic range, low background, and long lifetime
- DIFFRAC^{plus} - complete suite of WINDOWS NT based software

D8 ADVANCE – DON'T GAMBLE WITH YOUR ANALYTICAL RESULTS!




**find out
what's inside**

Germany:
Tel. (+49) 7 21/595-28 88
Fax (+49) 7 21/595-45 87
USA:
Tel. (+1) 608/276-3000
Fax (+1) 608/276-3006

BRUKER ANALYTICAL X-RAY SYSTEMS



Quantitative XRD



SIROQUANT™ Quantitative XRD Phase Analysis Software

Based on the Rietveld Method

SIROQUANT™ is a user friendly quantitative XRD analysis software package for the Windows™ operating system. It is suitable for the analysis of data files from any automated diffractometer. This is a commercially supported package and is used in 25 countries around the world.

SIROQUANT™ was developed to enable quantitative analysis of crystalline phases in geological and manufactured materials. The software can be used for the analysis of amorphous content of the sample and crystallite size. The software includes a database of over 300 common phases. The software also incorporates the analysis of clays and poorly crystalline materials.

The program supports the refinement of the following Rietveld parameters for each phase:


- scale factor
- preferred orientation
- unit cell dimensions
- halfwidths (U, V, W) (with a second halfwidth if required)
- the Pearson "m" and Pseudo-Voigt lineshape parameter

Refinement of the instrument zero and the line asymmetry parameter for the whole pattern is also supported.

SIROQUANT™ provides full profile Rietveld XRD analysis of crystalline phases in:

**Cements, Minerals, Raw Coal and Coal Ash,
Metals, Clays, Beach and Industrial Sands and
Inorganic Materials.**

Developed in conjunction with CSIRO in Australia, SIROQUANT™ is available directly from Sietronics Pty Ltd or through your local Bruker AXS office.



Sietronics


Analytical X-ray Systems

Sietronics Pty Limited ACN 008 536 794
14 Wales Street
Belconnen, ACT 2617, Australia
PO Box 3066
Belconnen, ACT 2616, Australia
Telephone +61-2-6251 6611
Facsimile +61-2-6251 6659
E-Mail sietron@ibm.net
Internet <http://www.sietronics.com.au>



Windows™ is a trademark of Microsoft Corporation.
All trademarks are the properties of their respective owners.

No matter what
it's made of, an X'Pert
can recognize it



With the Philips X'Pert X-ray diffractometer line, you'll find a solution to even the most complex materials identification problems. X'Pert-MRD for advanced research, X'Pert MPD for powder diffraction, stress/texture, microstructure and crystallographic studies.

A comprehensive hardware/software range offers unique features for ultimate accuracy, sensitivity, flexibility and ease of use.

For more information:

Philips Analytical

Lelyweg 1
7602 EA Almelo
The Netherlands
Tel. : +31 (546) 534444
Fax : +31 (546) 534598

Philips Analytical

85 McKee Drive, Mahwah, NJ
07430-1921, USA
Tel. : +1 (201) 5296111
Fax : +1 (201) 5295084

Internet:

<http://www.analytical.philips.com>

Let's make things better.

Published online by Cambridge University Press



PHILIPS

Editor-in-Chief

Deane K. Smith
1652 Princeton Drive
State College, PA 16803-3257 U.S.A.

Associate Editor-in-Chief**Editor for New Diffraction Data**

Ting C. Huang
6584 Radko Drive
San Jose, CA 95119-1924 U.S.A.

Managing Editor

Mary M. Rossi
JCPDS-International Centre for Diffraction Data
12 Campus Blvd.
Newtown Square, PA 19073-3273 U.S.A.

Western European Editor

Norberto Masciocchi
Università di Milano
Dipartimento di Chimica
Strutturale e Stereochimica Inorganica
Via Venezian, 21
20133 Milano, Italy

Eastern European Editor

Jaroslav Fiala
SKODA Research Ltd.
31600 Plzeň, Czech Republic

Editor for Australia and New Zealand

Brian H. O'Connor
Curtin University
GPO Box U 1987, Perth 6001
Western Australia, Australia

Editor for Japan

Hideo Toraya
Ceramics Research Lab
Nagoya Institute of Technology
Asahigaoka, Tajimi 507 Japan

International Reports Editor

Elizabeth Rivette Bond
17300 Seventeenth St. #J-133
Tustin, CA 92780 U.S.A.

Editorial Advisory Board

P. Bayliss, *Sydney, Australia*
C. Z. Bojarski, *Katowice, Poland*
A. Brown, *Dorset, England*
D. Cox, *Upton, New York*
W. Eysel, *Heidelberg, Germany*
J. Fiala, *Plzeň, Czech Republic*
L. Frevel, *Midland, Michigan*
P. Gado, *Budapest, Hungary*
H. Goebel, *Munchen, Germany*
R. Jenkins, *Newtown Square, Pennsylvania*
G. G. Johnson Jr., *State College, Pennsylvania*
Q. Johnson, *Livermore, California*
J. I. Langford, *Birmingham, U.K.*
D. Louër, *Rennes, France*
G. J. McCarthy, *Fargo, North Dakota*
H. F. McMurdie, *Gaithersburg, Maryland*
M. E. Mrose, *Gaithersburg, Maryland*
M. Nichols, *Livermore, California*
R. L. Snyder, *Columbus, Ohio*
E. K. Vasil'ev, *Irkutsk, Russia*
T. Yamanaka, *Tokyo, Japan*
R. A. Young, *Atlanta, Georgia*

AIP Production: Lin Miller, *Editorial Supervisor*

Andrea Witt, *Journal Coordinator*

Kelly Quigley, *Senior Production Editor*

Submit manuscripts (3 copies) to the most appropriate *Powder Diffraction* Editor listed on this page. The Editors will consider all manuscripts received, but assume no responsibility regarding them. Materials will be returned only when accompanied by appropriate postage. There is no publication charge. See *Powder Diffraction Notes for Authors* for additional information.

Powder Diffraction is a quarterly journal published by the JCPDS-International Centre for Diffraction Data through the American Institute of Physics (AIP). *Powder Diffraction* is a journal of practical technique, publishing articles relating to the widest range of application—from mineral analysis to epitaxial growth of thin films and to the latest advances in software. Although practice will be emphasized, theory will not be neglected, especially as its discussion will relate to better understanding of technique.

Proofs and all correspondence concerning papers in the process of publication should be addressed to: Editorial Supervisor, *Powder Diffraction*, AIP, 500 Sunnyside Blvd., Woodbury, NY 11797-2999.

For advertising rates and schedules contact AIP Advertising Department. Orders, advertising copy, and offset negatives should be sent to: Advertising Department, American Institute of Physics, 500 Sunnyside Blvd., Woodbury, NY 11797-2999; phone: (516) 576-2440; fax: (516) 576-2481.

Subscription Prices (1998)

	U.S.A & Canada	Mexico, Central & South America	Europe, Mid-East & Africa*	Asia & Oceania*
Individual	\$60	\$85	\$85	\$85
Institutional or Library	\$105	\$105	\$105	\$105

*Subscription rates to Eastern Hemisphere include air freight service.

Back-Number Prices. 1998 single copies: \$30. Prior to 1998 single copies: \$30.

Subscription, renewals, and address changes should be addressed to *AIP Circulation and Fulfillment Division (CFD)*, 500 Sunnyside Blvd., Woodbury, NY 11797-2999. Allow at least six weeks advance notice. For address changes please send both old and new addresses and, if possible, include a mailing label from the wrapper of a recent issue.

Claims, Single Copy Replacement and Back Volumes: Missing issue requests will be honored only if received within six months of publication date (nine months for Australia and Asia). Single copies of a journal may be ordered and back volumes are available in print or microform. Individual subscribers please contact AIP Circulation and Fulfillment Division (CFD) at (516) 576-2288; (800) 344-6901. Institutional or library subscribers please contact AIP Subscriber Services at (516) 576-2270; (800) 344-6902.

Reprint Billing: Contact: AIP Circulation and Fulfillment Division, Woodbury, NY 11797-2999; (516) 576-2234; (800) 344-6909.

Copying: Single copies of individual articles may be made for private use or research. Authorization is given (as indicated by the Item Fee Code for this publication) to copy articles beyond the use permitted by Sections 107 and 108 of the U.S. Copyright Law, provided the copying fee of \$6 per copy per article is paid to the Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, USA. Persons desiring to photocopy materials for classroom use should contact the CCC Academic Permissions Service. The Item Fee Code for this publication is 0885-7156/96 \$6.00.

Authorization does not extend to systematic or multiple reproduction, to copying for promotional purposes, to electronic storage or distribution, or to republication in any form. In all such cases, specific written permission from AIP must be obtained.

Permission for Other Use: Permission is granted to quote from the journal with the customary acknowledgment of the source. To reprint a figure, table, or other excerpt requires the consent of one of the authors and notification to AIP.

Requests for Permission: Address requests to AIP Office of Rights and Permissions, 500 Sunnyside Boulevard, Woodbury, NY 11797-2999; Fax: 516-576-2499; Telephone: 516-576-2268; E-mail: rights@aip.org.

Document Delivery: Copies of journal articles can be ordered for online delivery from the new Articles in Physics online document delivery service (URL: <http://www.aip.org/articles/>).

Reprints: Reprints can be ordered with or without covers only in multiples of 50 (with a minimum of 100 in each category) from AIP, Circulation & Fulfillment/Reprints, 500 Sunnyside Boulevard, Woodbury, NY 11797-2999; Fax: 516-349-9704; Telephone: 800-344-6909 (in U.S. and Canada), or 516-576-2234.



Powder Diffraction (ISSN: 0885-7156) is published quarterly (4X annually) by the JCPDS-International Centre for Diffraction Data through the American Institute of Physics. 1998 Subscription rates: US\$105. POSTMASTER: Send address changes to *Powder Diffraction*, AIP Circulation & Fulfillment Division, 500 Sunnyside Blvd., Woodbury, NY 11797-2999. Periodicals postage paid at Woodbury, NY 11797, and additional mailing offices.

Online Availability: Abstracts of journal articles published by the AIP and Member Societies (and several other physics publishers) are available in the SPIN database via the AIP online service PINET. Also available on PINET are *Advance Abstracts*, a current awareness service, other physics information resources, as well as internet services (Internet: elecprod@aip.org).

Copyright © 1998 JCPDS-International Centre for Diffraction Data, 12 Campus Blvd., Newtown Square, PA 19073-3273. All rights reserved.

THE FASTEST & EASIEST WAY TO
ACCOMMODATE YOUR TEXTURE,
STRESS, THIN FILM,
& POWDER APPLICATIONS

RIGAKU TAKES YOU FROM BRAGG-BRENTANO
TO PARALLEL BEAM AND BACK WITH NO ALIGNMENT



D/MAX 2000 PC
X-RAY DIFFRACTOMETER

NO NEED FOR DUAL 2θ ARMS
NO NEED TO REPLACE WHOLE COMPONENTS



Rigaku

RIGAKU/USA, INC. • 199 Rosewood Drive • Danvers, MA 01923


Phone: 978-777-2446 • Fax: 978-777-3594 • Website: www.rigaku.com

Deane K. Smith	Editorial: Changing of the Guard	195
J. A. Henao, J. M. Delgado, and M. Quintero	X-ray powder diffraction data and structural study of Fe_2GeSe_4	196
J. A. Henao, J. M. Delgado, and M. Quintero	X-ray powder diffraction data and structural study of Cd_4GeSe_6	202
Matteo Leoni, Paolo Scardi, and J. Ian Langford	Characterization of standard reference materials for obtaining instrumental line profiles	210
Yoshikazu Suzuki, Peter E. D. Morgan, and Koichi Niihara	Use of a high X-ray flux instrument for a mineral: X-ray powder diffraction pattern of $\text{CaMg}(\text{CO}_3)_2$	216
A. C. Stergiou, D. K. Giourelis	Crystallographic study of the mineral ferroan clinozoisite from three localities of Northern Greece	222
D. Mazza, M. Lucco-Borlera, and S. Ronchetti	Powder diffraction study of arsenic-substituted nasicon structures $\text{MeZr}_2\text{As}_{(3-x)}\text{P}_x\text{O}_{12}$ ($\text{Me}=\text{Na}^+, \text{K}^+$)	227
W. Wong-Ng, J. A. Kaduk, and W. Greenwood	Crystal structures and reference X-ray powder diffraction patterns of $\text{Sr}_{4-x}\text{Ca}_x\text{Pb}_2\text{O}_8$ ($x=1,2,3$)	232
Jialin Yan, Shiwei Wu, Xiangli Ou, Lingmin Zeng, and Jianmin Hao	X-ray powder diffraction data and structure refinement of CeFeGe_3	241
A. M. Xuriguera, C. Colbeau Justin, A. Elfakir, and M. Quarton	Crystal and powder neutron diffraction data of $\text{Ti}_3\text{Li}(\text{MO}_4)_2$ ($\text{M}=\text{Mo}, \text{W}$)	244
Nubuo Ishizawa, Atsushi Saiki, Kyoji Ohdan, and Mamoru Ai	Synthesis and powder X-ray diffraction data of a new iron phosphate $\text{Fe}(\text{PO}_4) \cdot 0.5\text{H}_2\text{O}$	246
J. Collazo, H. Hermann, A. Teresiak, and K. Wetzig	The effect of particles on the absorption factor for diffracted X-rays in single-phase powders measured at cylindrical geometry— a computational study	249
	International Reports	254
	Regional Report	254
	Computer Comments	255
	Calendar of Meetings	259
	Short Courses & Workshops	262
	Announcement	263
	Index to Volume 13	
	Subject Classification Scheme Used in Volume 13	264
	Subject Index to Volume 13	265
	Author Index to Volume 13	269

HUBER

Power in Powder Diffraction!

IMAGE FOIL GUINIER CAMERA 670



Get your diffractogram data within a couple of minutes! And get them without those bulky $K\alpha_2$ peaks.

Just have a look at the scan above. This is a zoomed view of sample 1g of the Commission on Powder Diffraction - Round Robin 1998. The dotted line was scanned on a conventional Bragg-Brentano instrument. It took some 4 hours for 100° 2-Theta with 3 sec per 0.02° step. The unbroken line results from our G670 and took 30 minutes for the same scan range but with 0.005° step size. Not to mention that $K\alpha_2$ is obsolete.

The samples are outside of the closed, and thus evacuable, camera housing. The Guinier Camera 670 fits onto the Basis 601 and uses any of the Guinier Monochromators of the HUBER Series 611/15. There are options for the Flat Specimen Sample Oscillation Device 670.1, or for the Capillary Rotation Device 670.2 together with the Furnace 670.3 and the High Temperature Controller 9634. The system includes the data collection software running on a PC. Data files of various popular ASCII types are available for further handling.

HUBER

X-RAY DIFFRACTION EQUIPMENT

HUBER Diffraktionstechnik GmbH
Sommerstrasse 4
D-83253 Rimsting, Germany

Tel.: +49 (0) 80 51 - 68 78 - 0
Fax: +49 (0) 80 51 - 68 78 - 10
e-Mail: Info@XHuber.com

crystallographica searchmatch

The makers of *Crystallographica* are proud to announce the launch of *Crystallographica Search-Match*, an all-new search-match program for Windows 95/98/NT.

- ◆ Works with all versions of the Powder Diffraction File including the new cPDF
- ◆ Search using full powder diffraction pattern and/or peak list
- ◆ Automatic residual search for multi-phase identification
- ◆ Unique integrated Boolean card retrieval and display
- ◆ Single / multi-phase full pattern powder simulations
- ◆ Reads common file formats
- ◆ Built-in tools include peak finding and background / $k-\alpha_2$ stripping
- ◆ Report writing directly to Microsoft Word
- ◆ Full 32-bit technology delivering unrivalled speed and power

Visit our new web site for details or contact us for a free demo CD!



OxfordCryosystems

3 Blenheim Office Park, Lower Road, Long Hanborough
Oxford · OX8 8LN · UK

Tel: +44 (0)1993 883488 · Fax: +44 (0)1993 883988

E.mail: info@OxfordCryosystems.co.uk

Search-Match? Upgrade to the state-of-the-art

www.crystallographica.co.uk

Mini Materials Analyser

High Performance Powder Texture Stress Thin Film Analyser

**2 Axis Harmonic Goniometer, 0.002 or 0.00002 deg. min. step, Spinning Powder Stage
Slim-line Eulerian Cradle (Max 2 Theta 160 deg.), Residual Stress.
Glancing Incidence Attachment, Thin Film analysis, X-ray Fluorescence analysis.**







**Use with
any X-ray
Generator!**

**Light-weight
benchtop
mounting!**

**Total access
lift-up cover!**

Software for Powder, Texture, Residual Stress
Crystallography, Reflectivity, Thin Film
Energy Dispersive X-ray Fluorescence
Qualitative analysis, All for Windows.



Visual Texture

Comment: Texture Data Collection Software

Control file: demo .dls

Output File: demo .sav

Count Time: 10

Pole	H	K	L	Norm. A	Norm. B	R.F.
1	1	0	0	0.000000	0.000000	26.562
2	1	1	1	0.000000	0.000000	45.700
3	1	0	2	0.000000	0.000000	41.700
Final						26.000

Status:

2THETA: 20

CHI: 0

PHI: 160

Counts: 2950

Diffraction Technology


Diffraction Technology Pty. Ltd.
Contact : Rod Clapp
Phone: 61 2 6242 8233
FAX: 61 2 6242 8266
Internet: rod@pcug.org.au

PO Box 444, Mitchell ACT 2911
Canberra, Australia.

United States : Jerry Krause
Analytical X-ray Instrument Sales
Phone/FAX : 303 277 1569
Internet: jkxray@aol.com


United Kingdom : Arthur Bagley
Hiltonbrooks Ltd.

Phone : 0 1477 534140
Fax : 0 1477 535041
Internet: arthur@xrays.u-net.com



Diftech 122D Automation System , XRD, XRF, Texture

Windows is a registered trademark of the Microsoft Corporation. Other product names are registered trademarks of their respective owners.



Performs Great Under Pressure

SPEX CertiPrep introduces a unique manual press which *performs with the greatest of ease!* Try one out today and start reducing your sample pelletizing costs to mere **peanuts** a day!

The Air-actuated Bench-Press™ eliminates the tedious task of hand pumping, necessary with other manual presses. Just push the button and the Bench-Press will generate up to **25 Tons** of pressure for sample pelletizing.

SPEX CertiPrep also offers the X-Press Series of 35 ton, automated laboratory presses, the 3630 and the 3624B.

A variety of pellet dies, reinforcing caps (SPEC-Caps), and Prep-Aid™ (sample binder) are also available to meet all of your sample pelletizing needs.



SPEX
INC.
CertiPrep 

203 Norcross Avenue,
Metuchen, NJ 08840

Website: www.spexcsp.com
Phone: 732-549-7144
Fax: 732-603-9647

1-800-LAB-SPEX • sampleprep@spexcsp.com (e-mail)

Add FLUXability to

Your Sample Preparation Techniques!


SPEX CertiPrep-Claisse automated fusion fluxers, "First and Finest in Fusion", quickly prepare hard-to-dissolve samples for XRF or ICP/AA analysis. Cement, slag, glass, ceramics, rocks and ores are fused with lithium borate fluxes to make glass discs or solutions. Accuracy, reliability and low cost per fusion are obtainable with the 3-burner and 6-burner models.



800-LAB-SPEX
www.spexcsp.com

SPEX
INC.
CertiPrep 

203 Norcross Avenue
Metuchen, NJ 08840
732-549-7144 Fax 732-603-9647
sampleprep@spexcsp.com



Rapid particle size reduction for XRD, XRF and IR



MCCRONE MICRONISING MILL


- Reproducible mean particle size - narrow distribution
- Unique grinding action (not a ball mill)
- Wet grinding minimises damage to sample crystallinity
- Choice of grinding elements (agate or corundum)

McCrone Scientific Ltd.

McCrone House, 155A Leighton Rd.
London NW5 2RD. U.K.
tel: 0171 267 7199 fax: 0171 267 3383

USA Distributor:

McCrone Accessories & Components
850 Pasquinelli Drive
Westmont, Illinois 60559-5531
phone: 630 887 7100 fax: 630 887 7764
mac@mccrone.com



END & SIDE WINDOW X-RAY PROPORTIONAL COUNTERS

LND manufactures a wide range of End & Side Window X-Ray Proportional Counters for your scientific and instrumentation applications.

LND's dependable and operationally proven detectors are designed and engineered for a long and stable counting life. Our exacting manufacturing procedures and strict, audited Quality Assurance policies meet DCAS MIL-Q 985A, MIL-E-1, and Appendix B of 10CFR50 Quality Control Standards. All LND detectors pass a rigid, computerized testing procedure. It's your assurance of unequalled performance. LND is your source for the highest quality radiation detectors. Call or write for information.



LND, INC.

3230 Lawson Blvd. Oceanside, N.Y. 11572 U.S.A.
Tel: (516) 678-6141 • Fax: (516) 678-6704 • Telex: 14-4563

TOM BAUM SERVICE & SALES
6072 FM 1390
KAUFMAN, TEXAS 75142
972-452-3116 VOICE & FAX

RE: X-RAY DIFFRACTION

Periodic Maintenance

Service contracts for parts and labor on X-ray Systems
Periodic Alignment of Goniometers.
Installation and training


X-RAY Component Repair:

Rebuild of Goniometers (includes Imron painting)
Scintillation Detector X-Tal replacement
High Voltage Tank repairs. vacuum pumped & oil refilled

UPGRADES and MODIFICATIONS

Equipped Machine shop for special designed parts
PC Automation. for Data acquisition and Analysis of
Diffraction Scans . DOS & Mac Systems available
Safety Upgrades. new style auto shutters. beam tunnels.
and slit assemblies.
Sample Changers. Sample Rotation. Monochromaters.
X-ray Tubes. etc.


XRD XRF EDS WDS SEM



MPD Multi-Purpose X-ray
Diffractometer


Versatile diffractometers specifically designed to take advantage of the rapid, real time data collection our advanced CPS X-ray detectors offer.

- * No scanning feature - acquire up to $120^\circ 2\theta$ simultaneously
- Windows 95 based user friendly software
- Unique capillary devices for analysis of air sensitive materials
- Identify materials even if only micro amounts are available - quickly and accurately



Capillary Furnace

- Microcrystalline polymorph identification
- Newly synthesized materials identification
- Thin films
- Electrochemistry
- Kinetic transformations of ceramics
- In-Situ Measurements



3D Plot 2θ vs. Temperature


for information on our complete product line please contact us

inel - Z.A - CD 405 - 45410 ARTENAY (FRANCE)
Tel. (33) 0 2 38 80 45 45 Fax. (33) 0 2 38 80 08 14
E.MAIL: inel@valcofim.fr-INTERNET: <http://www.valcofim.fr/inel>

inel Inc. P.O. Box 147, STRATHAM, NH 03885 (USA)
TEL. (603) 778-9161 FAX. (603) 778-9171
E-MAIL: inelinc@aol.com

Your complete source for:

- * X-ray Tubes
- * Capillaries (Quartz & Glass)
- * Goniometer Heads
- * Kodak™ X-ray film products
- * Huber Components
- * Complete Diffractometers
- * Classical X-ray Cameras
- * Custom Designs



Huber 424/511.1 Sealed Source Beamline

More than 4500 professionals at over 1400 research locations have chosen Blake Industries for their X-ray diffraction instrumentation.

Let us put our experience to work for you.



BLAKE INDUSTRIES, INC.


Blake Industries, Inc.

660 Jerusalem Road, Scotch Plains, New Jersey 07076 USA

Tel (908) 233-7240 - Email: blake4xray@worldnet.att.net - Fax (908) 233-1354

Your Software For X-Ray Analysis

Now optimized for the new, enhanced PDF database with more than 37,000 additional calculated patterns



Software package for data acquisition and analysis:

Diffractometer control, plot modules, peak search and calibration, qualitative and quantitative phase analysis, profile fitting and pattern simulation, lattice parameter refinement, crystallite size and micro stress.

Ready to plug to your PHILIPS and SIEMENS diffractometer hardware.



Stand alone package for qualitative and quantitative phase analysis.

Extremely fast search algorithms for ICDD-PDF database and user created databases. Very exact quantitative analysis using the WHOLE PATTERN method.




Extension package for most third party X-ray software. Graphic-analytical, profile-analytical and crystallographic options.

a. wassermann

röntgenanalytik.messsysteme.software

P.O. BOX 2631, D-87416 Kempten, Germany, Tel. (49)0831-79515, FAX (49)0831-79930, E-mail: RMSKempten@aol.com



**Improving
your IQ
has never
been easier !**

Identification & Quantification

EXPLORE THE NEW KRATOS/SHIMADZU X-RAY LABORATORY
SUPERIOR PERFORMANCE & PRODUCTIVITY
FOR THE IDENTIFICATION & QUANTIFICATION OF MATERIALS



XRF-1700

SEQUENTIAL X-RAY FLUORESCENCE
SPECTROMETER

**THE ONLY IMAGING
WAVELENGTH DISPERSIVE
INSTRUMENT IN THE WORLD..**

A Thin Window
4kW X-ray Tube Delivers
Unrivalled Power and
Performance For The
Detection Of Elements Be to U

Ask about our new low cost, benchtop EDXRF system



XRD-6000

PRECISION ENGINEERED X-RAY
DIFFRACTION SYSTEM

**THE MOST AFFORDABLE
MULTI-FUNCTIONAL
XRD INSTRUMENT AVAILABLE**

Designed For
Fast & Accurate
Identification & Quantification
Of Crystalline Compounds
And Structures

*Ask about our new applications booklet on
solid state pharmaceutical products*



Call or write today for more information:

PHONE: (914)426-6700 FAX: (914)426-6192

E-MAIL: XRD SYSTEMS: xrd@kratos.com

XRF SYSTEMS: xrf@kratos.com

INTERNET: <http://www.kratos.com>

KRATOS
ANALYTICAL
A SHIMADZU GROUP COMPANY

Expand your x-ray diffraction capabilities — not your overhead.

Whether you just need services or your XRD lab is on "overload," IC Laboratories provides every testing service and advanced capability you need in qualitative or quantitative x-ray diffraction analysis — from austenite to zeolites, from air filters to thin films. You are assured of rapid turn-around of results — as little as 48 hours — because IC Labs is one of the most highly automated commercial labs in the U.S., with knowledgeable personnel ready to address all your applications. For a copy of our technical prospectus, contact IC Laboratories.

IC Laboratories

Post Office Box 721
Amawalk, New York 10501
(914) 962-2477 www.ICLABS.com

We're the Specialists in XRD

Get specialty beryllium prototypes with production processing in mind.

Driven to perform? Come to us. We'll put the metal to the metal.

When your production schedule is short, you want to get it done right the first time. You can eliminate a major uncertainty when you entrust your custom beryllium fabrication, joining and coating to us. We'll help you shine.

Here's what we can do for you:

- Depend on us as the only fully integrated source for beryllium sheet and foil products. All critical operations are performed in-house for seamless excellence.
- We're your partner in joining and coating metals. When your project involves joining a metal to a metal, you can count on us.
- Come to us for all your UHV beryllium product needs: x-ray windows, chambers, beam pipes and more.
- Rely on our expert engineering. We like challenges. Taking your design concept into reality would be our pleasure.

Count on us for your tough jobs. Keep in mind that we're the only fully integrated source for beryllium products for the analytical, medical, and scientific industries. Challenge us with your "cutting edge" requirements. We look forward to hearing from you.

TEL: 510-623-1500 • FAX: 510-623-7600

E-Mail: Electrofusion@BrushWellman.com

44036 South Grimmer Boulevard • Fremont, California 94538 • USA



BRUSHWELLMAN

ELECTROFUSION PRODUCTS

We're the beryllium window folks.

Phase Analysis

Reflectometry

Texture Analysis


Stress Analysis

Max-Flux[®] Optics for X-ray Diffraction

High-Resolution Diffraction

Point-Defect Analysis

Thin-Film Analysis




Call today or visit our website to learn more about this innovative technology!

1788 Northwood Drive, Troy, Michigan 48084, USA
(248) 362-1290 • Fax (248) 362-4043 • <http://www.osmic.com>



Focusing lens with large collection angle and focal spot of down to 50 μm FWHM, accurately directs x-ray energies from 200 eV to 30 KeV.



Collimating lens collects large solid angle and produces low-divergence, well collimated beam.

New X-Ray Lenses For Rapid & Precise Measurements.

Need speed and accuracy in critical x-ray measurements? Specify the only commercially available polycapillary optics. Get real bottom line savings in time and money.

X-Ray Focusing

- ▲ Increased intensity
- ▲ Improved detection limits
- ▲ Greater spatial resolution

X-Ray Collimating

- ▲ Two dimensional collimation
- ▲ More efficient use of source for increased total flux
- ▲ Broad spectral bandwidth

Diffraction Applications

- ▲ Texture
- ▲ Stress
- ▲ Phase identification


Fluorescence Applications

- ▲ Composition
- ▲ Film thickness

Have a polycapillary optics application? Call Phil Bly, our Sales and Applications Engineer. It could be the most productive call you'll make this year.



30 Corporate Circle • Albany, New York 12203 • Phone: 518-464-3334 • FAX: 518-464-3335
info@xos.com www.xos.com



Publishing Agreement

1. It is our custom to request authors to vest the worldwide copyright of their papers with the JCPDS—International Centre for Diffraction Data for the full term of copyright and we would be grateful if you would confirm your acceptance of these terms by signing and returning the agreement below. We will not withhold permission for any reasonable request from you to publish any part of this paper in connection with any other work by you, provided the usual acknowledgements are given regarding copyright notice and reference to the original publication.
2. If it is appropriate, the author's employer may sign this agreement, and in any event the employer may reserve the right to use the paper internally or for promotional purposes only, by so indicating on this agreement. It is understood that proprietary rights other than copyright (including patent rights) are reserved.
3. If the author is a U. S. Government employee and this work was done in that capacity the assignment applies only to the extent allowable by U. S. law. If at least one co-author is *not* a government employee, said author should sign the agreement.
4. If the author is an employee of the British Government then HMSO will grant a non-exclusive license to publish this paper in the Journal in any form or media provided British crown copyright and user rights (including patent rights) are reserved.
5. The author warrants that the manuscript is the author's original work, and has not been published before. (If excerpts from copyrighted works are included, the author will obtain written permission from the copyright owners and show credit to the sources in the manuscript.) The author also warrants that the article contains no libelous or unlawful statements, and does not infringe on the rights of others.
6. If the work was prepared jointly the author agrees to inform co-authors of the terms of the agreement, and to sign on their behalf.

With the qualifications listed above I assign to JCPDS—International Centre for Diffraction Data, the copyright of my paper entitled

for publication in *Powder Diffraction*.

_____ Signed

_____ Date

If the paper is rejected, this assignment is null and void.