

Volume 24, Number 1 February 2018

# Microscopy AND Microanalysis



CAMBRIDGE  
UNIVERSITY PRESS

ISSN 1431-9276

# Incomparable...

...and still innovating



## **DIATOME** diamond knives

ultra 45° • cryo • histo  
ultra 35° • histo jumbo  
cryo immuno • ultra sonic  
ultra AFM & cryo AFM  
trimtool 20 • trimtool 45  
trimtool 90

Over 40 years of development,  
manufacturing, and customer service

## **DIATOME U.S.**

P.O. Box 550 • 1560 Industry Rd. • Hatfield, Pa 19440  
Tel: (215) 412-8390 • Fax: (215) 412-8450  
email: [sgkcck@aol.com](mailto:sgkcck@aol.com) • [www.emsdiasum.com](http://www.emsdiasum.com)

### Free customer service


Sectioning tests with biological and material research specimens of all kinds. We send you the sections along with the surfaced sample, a report on the results obtained and a recommendation of a suitable knife. Complete discretion when working with proprietary samples.

### Re-sharpening and reworking service

A re-sharpened Diatome diamond knife demonstrates the same high quality as a new knife. Even knives purchased in previous years can continue to be re-sharpened. The knives can be reworked into another type of knife for no extra charge, e.g. ultra to cryo or 45° to 35°.

### Exchange service

Whenever you exchange a knife we offer you a new DiATOME knife at an advantageous price.



# Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

THE OFFICIAL JOURNAL OF

MICROSCOPY SOCIETY OF AMERICA  
MICROANALYSIS SOCIETY  
MICROSCOPICAL SOCIETY OF CANADA /  
SOCIÉTÉ DE MICROSCOPIE DU CANADA  
MEXICAN MICROSCOPY SOCIETY  
BRAZILIAN SOCIETY FOR MICROSCOPY AND MICROANALYSIS  
VENEZUELAN SOCIETY OF ELECTRON MICROSCOPY  
EUROPEAN MICROBEAM ANALYSIS SOCIETY  
AUSTRALIAN MICROSCOPY AND MICROANALYSIS SOCIETY  
PORTUGUESE SOCIETY FOR MICROSCOPY

PUBLISHED IN AFFILIATION WITH

ROYAL MICROSCOPICAL SOCIETY  
GERMAN SOCIETY FOR ELECTRON MICROSCOPY  
BELGIAN SOCIETY FOR MICROSCOPY  
MICROSCOPY SOCIETY OF SOUTHERN AFRICA

## Editor-in-Chief

John Mansfield  
4304 Spring Lake Blvd.  
Ann Arbor, MI 48108-9657  
e-mail: thejfmjfm@me.com

## Administrative Editor

John Shields  
University of Georgia  
Athens, GA 30602  
e-mail: jpshield@uga.edu

## Biological Sciences Applications Editors

W. Gray (Jay) Jerome  
Department of Pathology, Microbiology and  
Immunology  
U-2206 MCN  
Nashville, TN 37232-2561  
e-mail: jay.jerome@vanderbilt.edu

Heide Schatten  
Department of Veterinary Pathobiology  
University of Missouri-Columbia  
Columbia, MO 65211  
e-mail: SchattenH@missouri.edu

Rosemary White  
CSIRO Plant Industry  
Canberra, ACT 2601, Australia  
e-mail: Rosemary.white@csiro.au

Elizabeth Wright  
Department of Pediatrics  
School of Medicine  
Emory University  
Atlanta, GA 30322  
e-mail: erwright@emory.edu

## Cultural Heritage Applications Editor

Edward P. Vicenzi  
Smithsonian Institution, Museum  
Conservation Institute  
4210 Silver Hill Rd., Suitland, MD 20746  
e-mail: VicenziE@si.edu

## Materials Sciences Applications Editors

Vinayak Dravid  
Materials Science and Engineering  
Northwestern University  
Evanston, IL 60208-3105  
e-mail: v-dravid@northwestern.edu

Georg E. Fantner  
Interfaculty Institute for Bioengineering  
École Polytechnique Fédéral de Lausanne  
Lausanne 1015, Switzerland  
e-mail: georg.fantner@epfl.ch

David J. Larson  
CAMECA  
5500 Nobel Drive  
Madison, WI 53711  
e-mail: david.larson@ametec.com

Ian MacLaren  
Materials and Condensed Matter Physics  
School of Physics and Astronomy  
University of Glasgow  
Glasgow G12 8QQ, UK  
e-mail: Ian.MacLaren@glasgow.ac.uk

Ross Marceau  
Institute for Frontier Materials  
Deakin University  
Geelong, VIC 3216, Australia  
e-mail: r.marceau@deakin.edu.au

Joseph Michael  
Sandia National Laboratories  
P.O. Box 5800, Albuquerque, NM 87185  
e-mail: jrmicha@sandia.gov

Yoosuf N. Picard  
Materials Science & Engineering  
Carnegie Mellon University  
Pittsburgh, PA 15213  
e-mail: ypicard@cmu.edu

Masashi Watanabe  
Dept. of Mater. Sci. & Eng.  
Lehigh University, Bethlehem, PA 18015  
e-mail: masashi.watanabe@lehigh.edu

## Special Issues and Reviews Editor

David J. Smith  
Department of Physics  
Arizona State University, Tempe, AZ 85287-1504  
e-mail: david.smith@asu.edu

## Book Review Editor

Cynthia Goldsmith  
Centers for Disease Control, Atlanta, GA 30333  
e-mail: csg1@cdc.gov

## M&M Program Guide Editor

Richard L. Martens  
1013 Bevell Building  
Box 870164, Tuscaloosa, AL 35487-0164  
e-mail: rmartens@caf.ua.edu

## Proceedings Editor

Gail Celio  
University of Minnesota, St. Paul, MN 55108  
e-mail: celio001@umn.edu

## Editorial Board

Ralph Albrecht	<i>University of Wisconsin, Madison, Wisconsin</i>
Ilke Arslan	<i>Pacific Northwest Laboratory, Richland, Washington</i>
Mary Grace Burke	<i>University of Manchester, Manchester, UK</i>
Barry Carter	<i>University of Connecticut, Storrs, Connecticut</i>
Wah Chiu	<i>Baylor College of Medicine, Houston, Texas</i>
Marc De Graef	<i>Carnegie Mellon University, Pittsburgh, Pennsylvania</i>
Niels de Jonge	<i>INM Institute for New Materials, Saarbrücken, Germany</i>
Elizabeth Dickey	<i>North Carolina State University, Raleigh</i>
Mark Ellisman	<i>University of California at San Diego, San Diego, California</i>
Pratibha Gai	<i>University of York, United Kingdom</i>
Marija Gajdardziska-Josifovska	<i>University of Wisconsin-Milwaukee, Milwaukee, Wisconsin</i>
Paul Kotula	<i>Sandia National Labs, Albuquerque, New Mexico</i>
William Landis	<i>University of Akron, Akron, Ohio</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Dale Newbury	<i>National Institute of Standards and Technology, Gaithersburg, Maryland</i>
Robert Price	<i>University of South Carolina, Columbia, South Carolina</i>
Jean-Paul Revel	<i>California Institute of Technology, Pasadena, California</i>
David Smith	<i>Arizona State University, Tempe, Arizona</i>
Nan Yao	<i>Princeton University, Princeton, New Jersey</i>
Nestor Zaluzec	<i>Argonne National Laboratory, Argonne, Illinois</i>

## Editorial Board Representatives from Affiliated Societies

Masashi Watanabe	<i>Lehigh University (MAS)</i>
Gautam Kumar Dey	<i>Bhabha Atomic Research Centre (EMSI)</i>
Gema Gonzalez	<i>Venezuelan Institute for Scientific Investigation (Venezuela)</i>
Michael Robertson	<i>Acadia University, Wolfville, Nova Scotia (Canada)</i>
Brendan Griffin	<i>University of Western Australia (AMMS)</i>
Guillermo Solorzano	<i>Pontificia Universidade Catolica, Rio de Janeiro (Brazil)</i>
Mike Matthews	<i>Atomic Weapons Establishment, Reading, Great Britain (EMAS)</i>
Miguel Yacaman	<i>Mexico Institute for Nuclear Research (Mexico)</i>
Henrique Almeida	<i>Universidade do Porto (Portugal)</i>

## Founding Editor

Jean-Paul Revel *California Institute of Technology, Pasadena, California*

## Previous Editors-in-Chief

Dale Johnson	<i>University of South Florida, Tampa, Florida</i>
Charles Lyman	<i>Lehigh University, Bethlehem, Pennsylvania</i>
Robert L. Price	<i>University of South Carolina, Columbia, South Carolina</i>

This journal is part of the **Cambridge Core** service. Access to online tables of contents and article abstracts is available to all researchers at no cost. Access to full-text articles online is provided to those with online subscription. Online subscriptions must be activated. Once your subscription is activated, free access to past, present, and forthcoming articles is available at:

***Microscopy and Microanalysis* website: [cambridge.org/MAM](http://cambridge.org/MAM).**


Instructions for authors submitting manuscripts may be found at [cambridge.org/MAM](http://cambridge.org/MAM). Select "Further Information" then select "Instructions for Contributors."

# TESCAN S8000G and novel Orage™ FIB column


*Anticipating the future, expanding your possibilities today*

## Orage™ FIB column: Prepare for the new level in sample preparation


- Next generation of Ga source FIB column and a guarantee of world-class quality in sample preparation
- Cutting-edge ion beam optics delivering improved resolution at low beam energies for damage-free ultra-thin TEM specimens
- High ion beam currents up to 100 nA enabling fast sputtering rates for maximum throughput and minimum time-to-result



Cross-section 50 µm wide prepared with an ion beam current of 85 nA.



1 keV FIB image of a lamella on the TEM grid.



50 µm-long cross-section prepared in a Li-ion battery cathode



S 8000

## Microscopy AND Microanalysis

*Microscopy and Microanalysis* publishes original research papers dealing with a broad range of topics in microscopy and microanalysis. These include articles describing new techniques or instrumentation and their applications, as well as papers in which established methods of microscopy or microanalysis are applied to important problems in the fields of biology or materials science. Microscopy and microanalysis are defined here in a broad sense, and include all current and developing approaches to the imaging and analysis of microstructure. The criteria for acceptance of manuscripts are the originality and significance of the research, the quality of the microscopy or microanalysis involved, and the interest for our readership.

Four types of communications are published in the Journal. **Regular Articles** are of substantial length and describe the findings of an original research project that satisfies the aims and scope of the Journal, described above. **Review Articles** summarize the current status of an important area within the aims and scope of the Journal. **Letters to the Editor** usually contain comments on recent articles that have appeared in the Journal. **Book Reviews** are also published, but these are solicited only through the Book Review Editor.

### Instructions for Contributors

Instructions for authors contributing manuscripts may be found at <http://mc.manuscriptcentral.com/mam> under "Resources: Instructions and Forms." Authors may also visit [cambridge.org/mam](http://cambridge.org/mam), select "Information," and then select "Instructions for Contributors."

### Copyright Information

Submission of a manuscript implies: that the work described has not been published before (except in the form of an abstract or as part of a published lecture, review, or thesis); that it is not under consideration for publication elsewhere; that its publication has been approved by all coauthors, if any, as well as by the responsible authorities at the institute where the work has been carried out; that, if and when the manuscript is accepted for publication, the authors agree to automatic transfer of the copyright to the Microscopy Society of America; that the manuscript will not be published elsewhere in any language without the consent of the copyright holders; and that written permission of the copyright holder is obtained by the authors for material used from other copyrighted sources.

All articles published in this journal are protected by copyright, which covers the exclusive rights to reproduce and distribute the article (e.g., as offprints), as well as all translation rights. No material published in this journal may be reproduced photographically or stored on microfilm, in electronic data bases, video disks, etc., without first obtaining written permission from the publisher.

The use of general descriptive names, trade names, trademarks, etc., in this publication, even if not specifically identified, does not imply that these names lack protection by the relevant laws and regulation.

Authorization to photocopy items for internal or personal use, or the internal or personal use of specific clients, is granted by Cambridge University Press, provided that the appropriate fee is paid directly to Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, USA [Tel: (508) 750-8400], stating the ISSN (1431-9276), the volume, and the first and last page numbers of each article copied. The copyright owner's consent does not include copying for general distribution, promotion, new works, or resale. In these cases, specific written permission must first be obtained from the publisher.

### Disclaimer

The Microscopy Society of America, the other societies stated, and Cambridge University Press cannot be held responsible for errors or for any consequences arising from the use of the information contained in this journal. The appearance of scientific reports and/or workshops, or any other material in *Microscopy and Microanalysis* does not constitute an endorsement or approval by The Microscopy Society of America of the findings, data, conclusions, recommendations, procedures, results, or any other aspect of the content of such articles. The appearance of advertising in *Microscopy and Microanalysis* does not constitute an endorsement or approval by

## Aims and Scope

The Microscopy Society of America of the quality or value of the products advertised or any of the claims, data, conclusions, recommendations, procedures, results, or any other information included in the advertisements.

While the advice and information in this journal is believed to be true and accurate at the date of its going to press, neither the authors, the editors, nor the publisher can accept any legal responsibility for any errors or omissions that may be made.

### Subscription Information

*Microscopy and Microanalysis* is published bimonthly in February, April, June, August, October, and December by Cambridge University Press. Three supplements (*Meeting Guide*, *Program Guide*, and *Proceedings*) are published in June and August.

**Society Rates:** Members of the Microscopy Society of America should contact the MSA Business Office for all subscription inquiries: Microscopy Society of America, 11130 Sunrise Valley Dr, Suite 350, Reston, VA 20191, Tel.: (703) 234-4115, Email: [associationmanagement@microscopy.org](mailto:associationmanagement@microscopy.org), URL: [www.microscopy.org](http://www.microscopy.org). Members of other affiliated societies should contact their respective society business offices for all subscription inquiries.

**Subscription Rates:** Institutions print and electronic: US \$2364.00 in the USA, Canada, and Mexico; UK £1422.00 + VAT elsewhere. Institutions online only: US \$1533.00 in the USA, Canada, and Mexico; UK £928.00 + VAT elsewhere. Individuals print plus online: US \$691.00 in the USA, Canada, and Mexico; UK £420.00 + VAT elsewhere. Prices include postage and insurance.

**USA, Canada, and Mexico:** Subscribers in the USA, Canada, and Mexico should send their orders, with payment in US dollars or the equivalent value in Canadian dollars, to: Cambridge University Press, Customer Services Department (Journals), 1 Liberty Plaza, New York, NY 10006, USA. Tel: (845) 353-7500. Fax: (845) 353-4141. Orders may be phoned direct (toll free): (800) 872-7423. E-mail: [journals\\_subscriptions@cup.org](mailto:journals_subscriptions@cup.org).

**Outside North America:** Subscribers elsewhere should send their orders, with payment in sterling, to: Customer Services Department (Journals), Cambridge University Press, University Printing House, Shaftesbury Road, Cambridge CB2 8BS, UK. Tel: +44 (0)1223 326070. Fax: 44 (0)1223 325150. E-mail: [journals@cambridge.org](mailto:journals@cambridge.org)

**Change of address:** Allow six weeks for all changes to become effective. All communications should include both old and new addresses (with postal codes) and should be accompanied by a mailing label from a recent issue. Society members should contact their respective society business offices to inform them of address changes.

### Editorial Office

John Mansfield, Editor in Chief, 4304 Spring Lake Blvd., Ann Arbor, MI 48108-9657, USA; Tel: (734) 994-3096; Fax: (734) 763-2282; E-mail: [thejfmjfm@me.com](mailto:thejfmjfm@me.com).

### Office of Publication

Cambridge University Press, 1 Liberty Plaza, New York, NY 10006, USA; Tel: (212) 337-5000; Fax: (212) 337-5959.

### Advertising Sales & Production

Kelly Miller, M.J. Mrvica Associates, Inc., 2 West Taunton Avenue, Berlin, NJ 08009, USA; Tel: (856) 768-9360; Fax: (856) 753-0064.

© 2018 by Microscopy Society of America. Printed in the United States on acid-free paper. Periodicals postage paid at New York, NY, and additional mailing offices. Return postage guaranteed. Postmaster: Send address changes in the U.S.A. and Canada to *Microscopy and Microanalysis*, Subscription Department, Cambridge University Press, 1 Liberty Plaza, New York, NY 10006.

# SEM & TEM Sample Preparation

## Microscopy



### DESK V TSC

Sample Coater

- Fully automated
- Minimize sample charge
- Reduce e-beam damage
- Increase secondary electron emissions

### BENCH TOP TURBO

Evaporation Tool for Tissue and Life Science

- Carbon coating
- Sputtering
- AC glow discharge
- Low-voltage thermal evaporation

## Support Performance Value

### DV-502

Large Sample, High Throughput Sputter & Evaporation

- Rotary shadowing for metals and carbon
- High vacuum gate valve for rapid cycling
- Recipe development

**DENTON VACUUM**

Enabling Innovation

All Denton Vacuum Microscopy tools are backed by a 12 month warranty, over 50 years of process knowledge, an in-house process engineering group, worldwide support, and Global factory service center.

To learn more about how to better prepare your samples, contact your electron microscope manufacturer or call us at +1-800-666-6004


[www.dentonvacuum.com](http://www.dentonvacuum.com)

# Nano-Positioning for Microscopy & Imaging



PIFOC® nano-focusing drive

**Piezo stages & positioners** are essential tools for high-resolution microscopy, such as Super Resolution Microscopy or AFM. Their sub-atomic resolution and extremely fast response allow researchers to create higher-quality images faster. PI provides a large variety of fast positioning stages and piezo objective nano-positioners for 3D imaging (Z-stack acquisition), fast-focusing applications, and light sheet microscopy.



Piezo flexure actuators excel in applications from fast nano-dosing to precision alignment.

Affordable XY & XYZ Piezo Stages for SR Microscopy: P-545 PInano®




**PI**

Physik Instrumente  
[www.pi-usa.us](http://www.pi-usa.us)  
 508-832-3456 (East)  
 949-679-9191 (West)



PRECISION | SPEED | STABILITY - MOTION CONTROL & POSITIONING SOLUTIONS



Thermo Scientific™ Avizo™ 3D visualization of two large adjacent crystalline dendrites of a bulk-metallic-glass matrix composite ( $Zr_{58.5}Ti_{14.3}Nb_{5.2}Cu_{6.1}Ni_{4.9}Be_{11.0}$ ). Data was obtained by large volume serial sectioning tomography using the Thermo Scientific™ Helios™ PFIB DualBeam™ microscope. The sectioned block is about  $90 \times 80 \times 70 \mu m^3$ . Sample from The University of Tennessee, USA. Images courtesy of The University of Manchester.

## Large 3D volumes with unprecedented surface resolution

Until recently, available technologies have limited the volumes and depths of materials that can be analyzed at high resolution, ultimately restricting the insight into structural, crystallographic, and chemical properties. This is no longer the case. The Thermo Scientific Helios PFIB DualBeam microscope offers unrivaled access to regions of interest deep below the surface—combining serial section tomography with statistically relevant data analysis.

[thermofisher.com/EM-materials-science](https://thermofisher.com/EM-materials-science)

**ThermoFisher**  
SCIENTIFIC



# Microscopy AND Microanalysis

An International Journal for the Biological and Physical Sciences

Volume 24, Number 1  
February 2018

## MATERIALS SCIENCE APPLICATIONS


- Quantification of Olivine Using Fe  $L\alpha$  in Electron Probe Microanalysis (EPMA)  
*Ben Buse and Stuart Kearns* 1
- The Influence of Beam Broadening on the Spatial Resolution of Annular Dark Field Scanning Transmission Electron Microscopy  
*Niels de Jonge, Andreas Verch, and Hendrix Demers* 8

## BIOLOGICAL SCIENCE APPLICATIONS

- A Simple Method for the Cross-Section Area Determination of Single Profiled Fibers and Its Application  
*Zhangchuan Peng, Chun Liu, Lin Zhang, Wei Li, Wenbo Hu, Sanyuan Ma, and Qingyou Xia* 17
- Overexpression of MicroRNA-29b Decreases Expression of DNA Methyltransferases and Improves Quality of the Blastocysts Derived from Somatic Cell Nuclear Transfer in Cattle  
*Shuang Liang, Zheng-Wen Nie, Jing Guo, Ying-Jie Niu, Kyung-Tae Shin, Sun A. Ock, and Xiang-Shun Cui* 29
- Early Effects of Ionizing Radiation on the Collagen Hierarchical Structure of Bladder and Rectum Visualized by Atomic Force Microscopy  
*Svetlana L. Kotova, Peter S. Timashev, Galina V. Belkova, Marina V. Kochueva, Ksenia V. Babak, Victoria A. Timofeeva, Elena B. Kiseleva, Olga O. Vasilieva, Anna V. Maslennikova, and Anna B. Solovieva* 38
- Superficial and Inner Examination of a Microwave-Irradiated Dental Acrylic Resin and Its Metal-Polymer Interface  
*Marian C. Popescu, Bogdan I. Bitu, Vasilica Tucureanu, Dan Vasilache, Melania A. Banu, Andrei M. Avram, and Raluca A. Giurescu-Dumitrescu* 49

## MICROGRAPHIA

- An Open-Source Storage Solution for Cryo-Electron Microscopy Samples  
*Eveline Ultee, Fred Schenkel, Wen Yang, Susanne Brenzinger, Jamie S. Depelteau, and Ariane Briegel* 60
- Determining Direction of Axonal Flow in the Equine Ramus Communicans by Ultrastructural Examination of the Plantar Nerves 2 Months after Transecting the Ramus  
*Fakhri Al-Bagdadi, Jim Schumacher, Jessi Carter, Ferenc Tóth, and Robert W. Henry* 64
- Is the Binding Pattern of Zinc(II) Equal in Different Bryophyte Species?  
*Marko S. Sabovljević, Marieluise Weidinger, Aneta D. Sabovljević, Wolfram Adlassnig, and Ingeborg Lang* 69
- Polarization Microscopy and Infrared Microspectroscopy of Integument Coverings of Diapausing Larvae in Two Distantly Related Nonsocial Bees  
*Maria Luiza S. Mello, Benedicto de Campos Vidal, and Jerome G. Rozen, Jr.* 75



**On the Cover:** Differential interference contrast view of a wrinkled integument waxy cover from a *Hesperapis rhodocera* larva. For further details, see Mello et al., pp. 75–81.

**Microscopy and Microanalysis website:** <http://www.cambridge.org/MAM>

Indexed in Chemical Abstracts, Current Contents, BIOSIS, and MEDLINE (PubMed)



Special tweezer slots allow easy and secure removal of Cryo Grid Boxes.

# EMS has it!

## Cryo-EM Pucks

- 12 wells per puck for Cryo Grid Boxes
- Each puck has a unique alpha-numeric code for easy identification
- Indexed wells for sample tracking
- Holds round Cryo Grid Boxes
- Puck depth accommodates pin type lid style Cryo Grid Boxes
- When using Cryo Grid Boxes with flat-style lids, you can store up to 24 Cryo Grid Boxes per puck

Simplify how you store and transport your cryogenically cooled Cryo-EM specimen grids. Each puck holds up to 24 round Cryo Grid Boxes (96 grids). A custom serial number on each puck allows easy sample tracking. A 10-puck capacity storage cane and a 7-puck capacity shelved shipping cane are also available.


# Electron Microscopy Sciences

P.O. Box 550 • 1560 Industry Rd.  
Hatfield, Pa 19440  
Tel: (215) 412-8400  
Fax: (215) 412-8450  
email: [sgkck@aol.com](mailto:sgkck@aol.com)  
or [stacie@ems-secure.com](mailto:stacie@ems-secure.com)

[www.emsdiasum.com](http://www.emsdiasum.com)

Organized storage and transport for Cryo-EM specimen grids under cryogenic conditions

# NEW: EMS Cryo Pucks




71168-01 — EMS Cryo Pucks Complete Set

(7) EMS Cryo Pucks, (1) Shelved Puck Shipping Cane, (1) Shelved Storage Cane, (1) Angled Cryo-Tongs, (1) Protective Storage Case

71168-09

Shelved Storage Cane



71168-10



71168-12



71168-13

## Ordering Information

Cat.-No.	Description	Qty.
71168-01	EMS Cryo Pucks Complete Set <b>Includes:</b> (7) EMS Cryo Pucks (1) Shelved Puck Shipping Cane (1) Shelved Storage Cane (1) Angled Cryo-Tongs (1) Protective Storage Case	set
71168-02	EMS Cryo Pucks	each
71168-07	EMS Cryo Pucks	7/pk

Cat.-No.	Description	Qty.
71168-08	EMS Cryo Puck Storage Case	each
71168-09	Shelved Storage Cane (holds 10 pucks)	each
71168-10	Double Puck Loading Dewar with Lid	each
71168-11	Bent Cryo Tong	each
71168-12	Shelved Puck Shipping Cane	each
71168-13	Cryo Express Dry Shipper with Case	each