

now available
from EMS...

the clever
cleaver
you **can't**
live without!


Meet the LatticeAx™ — the small,
accurate, fast, low-cost cleaving
solution, suitable for any lab.

its simple...


just set...



indent...



cleave



LatticeAx™ Cleaving System

The LatticeAx™ is a precision cleaver that fits in the palm of your hand. In an amazingly small footprint (4" cube, 100mm³), the patent pending "Ax" and process are designed to assist the user to cleave wafers, strips, or pieces to precisely sized samples with localized targets. Unlike automated cleaving systems, the LatticeAx™ bridges the gap between manual cleaving and expensive alternatives. With the smart combination of a high-magnification camera and precision knobs, the LatticeAx™ enables highly accurate, manual cleaving of wafer pieces and small samples of any aspect ratio.




Electron Microscopy Sciences

P.O. Box 550
1560 Industry Road
Hatfield, PA 19440
Tel: (215) 412-8400
Fax: (215) 412-8450
email: sgkcck@aol.com
or stacie@ems-secure.com
www.emsdiasum.com

applications

- Site specific or general area cross-section for SEM
- Target localization prior to FIB
- Prep for Broad Beam ion milling or polishing
- Downsize samples for SEM with height restriction
- Cleaving to create uniform samples for other analysis tools with non-wafer scale stages
- Vertical, mirror image cleaving for photonics analysis

Contact us about the full line of LatticeGear Wafer Cleaving Tools



Living up to Life

Leica

MICROSYSTEMS



Coat with Confidence

The **NEW** Leica EM ACE range of coaters
for EM sample preparation

It's the detail captured in the coater that becomes the detail captured in your image. Our **NEW Leica EM ACE Coaters** are designed to excel at just that. With innovations from design to interface, you can be sure that whatever element you coat with, you'll also be coating with the elements of **Simplicity**, **Reproducibility**, **Performance** and above all, **Confidence**.

www.leica-microsystems.com/mm-coaters



© 2013 Leica Microsystems, Inc. BGA#2013-2



SKYSCAN 1272: New generation of microtomographs from BRUKER MICROCT

EXCEPTIONAL PERFORMANCE BEYOND RESEARCHER'S IMAGINATIONS



free QR-reader app
is available in the AppStore

- Up to 209 Megapixel (14450 x 14450 pixels) in every virtual slice through objects
- More than 2600 such slices can be reconstructed after a single scan
- Due to phase-contrast enhancement, object details as small as 0.4 μ m can be detected
- Maximum scanning diameter - 75mm, integrated micropositioning stage
- Multithreaded CPU / GPU 3D reconstruction for single computer or cluster
- Possibility to scan samples during compression, tension, cooling or heating.
- 16-position sample changer with automatic adjustment of magnification and scanning protocol
- At the end of the scan the system can send you e-mail with a direct link to the scanning results
- The results can be exported to iPad or iPhone for 3D rendering by special software for mobiles

bruker-microct.com

microtomography

Innovation with Integrity

M&M 2013 MICROSCOPY & MICROANALYSIS

August 4 – 8 • Indianapolis, IN

Register
now—save
\$100!

SHARE

information and
techniques

LEARN

about cutting-edge
technology


EXPAND

your professional
network

GROW and

SOLIDIFY

your career!



Register Online: <http://microscopy.org/MandM/2013>

Visit M&M website often for detailed program information,
registration & hotel info.



International
Metallographic
Society
An Affiliate Society of ASM International®

Preparation Equipment and Microscopy Supplies

The single source for All your microscopy supplies and specimen preparation equipment.



- Vacuum Coating Systems
- Calibration Standards
- PELCO® easiGlow™ Glow Discharge Unit
- SEM Sample Holders and Mounts
- Silicon Nitride TEM Membranes
- PELCO BioWave Pro® Tissue Processor
- TEM Support Films



- AFM Supplies
- Quality Laboratory Tweezers
- Vacuum Pick-up Systems
- Digital Stereo Microscopes
- Conductive Adhesives
- FIB Supplies



Complete line of compact Cressington EM Sample Coaters.

 **TED PELLA, INC.**
Microscopy Products for Science and Industry


sales@tedpella.com

800-237-3526


www.tedpella.com

The moment your data change
scientific minds.
This is the moment we work for.

// RECOGNITION
MADE BY CARL ZEISS



Correlative microscopy with confocal and SEM of a butterfly wing. Green: reflection; red: auto fluorescence.



Correlative microscopy with superresolution dSTORM and SEM of yeast. Calcofluor (blue, cell wall); red: G protein coupled receptor (GPCR), localized in vesicular structures.

As a class-leading manufacturer of both light and electron microscopes, Carl Zeiss is uniquely positioned to enable correlative microscopy studies. Precisely overlay fluorescence images acquired with widefield, confocal or superresolution microscopy with the ultraresolution information collected with electron microscopy.


www.zeiss.com/corrmic



We make it visible.

Count on Us!


Best EDS Performance with the **NEW** Slim-line XFlash® 6



You can count on the **NEW** XFlash® SDD generation:


- Best solid angle – optimum geometry and active areas from 10 mm² to 100 mm²
- Best throughput – up to 600,000 cps output at 1,500,000 cps input
- Best energy resolution – 121 eV at Mn K α , 47 eV at F K α , 38 eV at C K α (FWHM, exceeds ISO 15632:2002 requirements)

www.bruker.com




PELCO BioWave® Pro

Fast, Low Temperature
Microwave Tissue Processing



- Rapid specimen turnaround
- Low temperature environment
- Programmable processing
- Consistent quality results



Application kits for paraffin processing, decalcification, confocal, light and electron microscopy available.



TED PELLA, INC.
Microscopy Products for Science and Industry

www.tedpella.com sales@tedpella.com 800.237.3526

Save the Date!

AUGUST 3–7, 2014 • M&M 2014 & IUMAS-6



held in conjunction with



The SPI Supplies Family of Instruments



Excellence in sample preparation
just a click away... 2spi.com

Your results will never be better than your sample preparation. See how SPI Supplies can help you deliver the highest quality results for all your SEM/EDS, TEM and FESEM applications.

1. **Osmium Plasma Coaters**
for FESEM Applications
2. **SPI-MODULE™**
Sputter/Carbon Coater Module
3. **Gentle Mill™**
Ion Milling System
4. **Plasma Prep™ III Solid State**
Plasma Cleaner for cleaning TEM holders

5. **Plasma Prep™ X**
Parallel Plate Plasma Etcher
6. **Plasma Prep™ III Plasma Etcher**
with PPIII Process Controller
7. **Vacu Prep™ II**
Turbo Pump Evaporation System
8. **SPI-DRY™**
Critical Point Dryer
9. **Precision Spin Coater**
Spin coater



SPI Supplies Division of **STRUCTURE PROBE, Inc.**

P.O. Box 656 • West Chester, PA 19381-0656 USA
Phone: 1-610-436-5400 • 1-800-2424-SPI (USA and Canada) • Fax: 1-610-436-5755 • 2spi.com • E-mail: sales@2spi.com



CAMBRIDGE


New to the Advances in Microscopy and Microanalysis book series!

Scanning Electron Microscopy for the Life Sciences

Heide Schatten

University of Missouri, Columbia

US\$120.00; Hb: 978-0-521-19599-7; 312 pp



Recent developments in scanning electron microscopy (SEM) have resulted in a wealth of new applications for cell and molecular biology, as well as related biological disciplines. It is now possible to analyze macromolecular complexes within their three-dimensional cellular microenvironment in near native states at high resolution, and to identify specific molecules and their structural and molecular interactions. New approaches include cryo-SEM applications and environmental SEM (ESEM), staining techniques and processing applications combining embedding and resin-extraction for imaging with high resolution SEM, and advances in immuno-labeling. With chapters written by experts, this guide gives an overview of SEM and sample processing for SEM, and highlights several advances in cell and molecular biology that greatly benefited from using conventional, cryo, immuno, and high-resolution SEM.

About the series

The Press currently publishes the Microscopy and Microanalysis (MAM) journal in conjunction with the MSA, which reaches 4,000 microscopists and is affiliated with 12 international microscopy societies. The series would be a natural development from this journal, and will take a broad view of the discipline, covering topics from instrumentation to imaging, methodology and analysis across physical science, materials science, biology and medicine. Books commissioned for the series will range from advanced undergraduate textbooks through to research and practitioner oriented monographs for researchers. The series aims to produce a coherent source of material, encouraging the communication and exchange of ideas across these divergent fields, ensuring that the series appeals to a broad community in the physical and life sciences.

Forthcoming titles in this series:

Microscopic Nanocharacterization of Materials

by Michael Isaacson

Energy Filtered Electron Microscopy and Electron Spectroscopy

by Richard Leapman

Dynamic and Transmission Electron Microscopy

by Nigel Browning, Thomas LaGrange, Bryan Reed, Henning Stahlberg, Bradley Siwick



www.cambridge.org/us

800.872.7423



CAMBRIDGE
UNIVERSITY PRESS

Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in optical, electron or scanning probe microscopy, the biological or the physical sciences, MSA takes your knowledge to the next level!

Members Receive:

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Plus discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and much more!








Join MSA Today!

For more information: visit www.microscopy.org or call 1-800-538-3672



High Definition Digital
TEM Cameras with
1 to 16 Megapixels

- AMT SOLUTIONS
- Life Science Cameras
- Material Science Cameras
- Easy To Use Software
- Reliability and Services
- TEM Integration
- Extensive Support



3 Electronics Avenue, Danvers, MA 01923 • Phone: 978.774.5550 • www.amtimaging.com

Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in light, electron, or scanning probe microscopy, or the biological or the physical sciences, MSA takes your knowledge to the next level!

Members Receive:




- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and more!



Join MSA Today!

Visit: www.microscopy.org or call 1-800-538-3672

**the highest quality...
the most precise sectioning...
incomparable durability**



Free Customer Service

Sectioning tests with biological and material research specimens of all kinds. We send you the sections along with the surfaced sample, a report on the results obtained and a recommendation of a suitable knife. Complete discretion when working with proprietary samples.

Re-sharpening and Reworking Service

A re-sharpened Diatome diamond knife demonstrates the same high quality as a new knife. Even knives purchased in previous years can continue to be re-sharpened. The knives can be reworked into another type of knife for no extra charge, e.g. ultra to cryo or 45° to 35°.

Exchange Service

Whenever you exchange a knife we offer you a new Diatome knife at an advantageous price.

40 years of development,
manufacturing, and
customer service

DiATOME
diamond knives







ultra 45° • cryo • histo • ultra 35° • STATIC LINE II
cryo-P • cryo immuno • ultra sonic
cryotrim 45 and 25 ultra • AFM & cryo AFM • cryo 25°

DiATOME US • P.O. Box 550 • 1560 Industry Rd. • Hatfield, Pa 19440

Tel: (215) 412-8390 • Fax: (215) 412-8450 • email: sgkcck@aol.com • www.emsdiasum.com

The Ion Mill You Know: Reborn

Low energy ion guns Cold stage Touch screen Automation and recipes
Live viewing and saving images in Gatan DigitalMicrograph®



Top images: Benefits of low voltage operation: ZnO on Sapphire. Images courtesy of P. Vernegués - CRHEA-CNRS, France.

Bottom images: Images at increasing magnification, all stored in Gatan DigitalMicrograph®, left from PIPS™ II optical camera, others from Gatan Orius™ SC200B camera on TEM, at 200 kV. SrRuO₃/5nmPZT/SrRuO₃ thin film on SrTiO₃ (001) with PIPS™ II final polishing step at 300 Volts.



www.gatan.com

PIPS II

Precision Ion Polishing System II

Enhancing Electron Microscopy


ANALYTICAL TEM DIGITAL IMAGING SPECIMEN PREPARATION

TEM SPECIMEN HOLDERS SEM PRODUCTS SOFTWARE

<https://doi.org/10.1017/S1431927613001797> Published online by Cambridge University Press

Ultra-versatility Sub-nm Resolution

JEOL introduces a new series of Field Emission Scanning Electron Microscopes with ultrahigh resolution and outstanding low kV/low vacuum performance for imaging and nanoanalysis of structures, surface details, and magnetic samples. From the budget-conscious lab to leading-edge research center, JEOL offers superior price/performance in the new 7000 series SEMs.



*See what you've been missing.
The new FE-SEM Series from
JEOL - available today.
jeolusa.com/FESEM*

JEOL

Another
Extreme Imaging
Solution

Global Solutions Provider for Advanced Technology
www.jeolusa.com • salesinfo@jeol.com
978-535-5900



7000 **FE-SEM** Series