


MICROSCOPY TODAY

MARCH 1995

ISSUE #95-2

<https://doi.org/10.1017/S1551929500063021> Published online by Cambridge University Press



MICROSCOPY TODAY

PO BOX 620122
Middleton, WI 53562

Bulk Rate
U.S. Postage
PAID
Permit No 1547
Madison, WI

Resolution 5nm at $E_{PE} = 500eV$

Resolution 4nm at $E_{PE} = 1keV$

Partners in Discovery:
Electron Microscopes from Carl Zeiss

Resolution 2.5nm at $E_{PE} = 5keV$

200 nm

Resolution 1.2nm at $E_{PE} = 20keV$

The advertisement features a central photograph of a woman in a patterned shirt standing next to a large, white scanning electron microscope (SEM). The background is a grid of SEM images showing various particles. Text boxes are overlaid on the images, providing resolution data for different accelerating voltages. A scale bar labeled '200 nm' is also present.

DSM 982 GEMINI The New-Technology SEM

If you are looking for a field-emission scanning electron microscope with this outstanding resolution, with fully quantitative EDX-analysis capability at high take-off angle and with the potential to tilt your 6 inch sample

to 50° at a working distance of only 6 mm, then you should consider our revolutionary DSM 982 GEMINI. In any field of application your results will far exceed anything ever attainable before.

Please call us for further information:

Carl Zeiss, Inc.
Electron Optics Division
Thornwood, NY 10594
Tel. (800) 356-1090
Fax (914) 681-7443



Carl Zeiss
Electron Optics Division
D-73446 Oberkochen
Germany
Tel.: (0 73 64) 20-27 00
Fax: (0 73 64) 20-45 30