







XEIA3


UHR SEM/i-FIB

First class scientific instrument



765 Commonwealth Drive, Suite 101
Warrendale, PA 15086
Tel: (724) 772 - 7433


tescan-usa.com




TESCAN
XEIA3




Visit the meeting website often for details on
Symposia, Exhibits, Biological and Physical Sciences Tutorials,
In-Week Intensive Workshops, Posters,
Sunday Short Courses and Pre-Meeting Congresses



Monocytes (White Blood Cells) Adhering to the Inside Surface of an Artery as Part of an Inflammatory Reaction. W. Gray (Jay) Jerome, Vanderbilt University



Cast A347 Alloy Made by Semi-solid Melting (Mert Fleming's Development) Weck's Reagent in Bright Field. George Vander Voort, Consultant (Struers Inc.)



High Density Lipoprotein (HDL; the good cholesterol carrier) Stacking Together in Solution. W. Gray (Jay) Jerome, Vanderbilt University



M&M 2016
MICROSCOPY &
MICROANALYSIS
July 24-28, 2016 • Columbus, Ohio

<http://microscopy.org/MandM/2016>

Program Information | Information for Students | Meeting Awards | Exhibitor List & Floor Plan
Registration (opens February 2016) | Columbus Hotel Information (reservations available February 2016)



Isn't it about time you had your own Digital Microscope?


Now you can with the portable, affordable uScopeMXII!

The uScopeMXII is a small digital desktop microscope you can use in your workplace or home office. It captures images from standard glass slides and sends them to your PC.

You can interactively browse slides with full control of focus, image processing, and location. You can also scan regions of interest creating fully focused image sets.

The industry-standard USB interface makes it simple to plug in and get started. It easily interfaces with your desktop or laptop PC and allows you to view and capture slide images in a wide variety of environments.


The uScopeMXII is designed and manufactured in the United States.



<http://www.uscopes.com/dz1>



<http://www.uscopes.com/dz2>



<http://www.uscopes.com/dz3>



**Call us today for a quote
or to schedule a demonstration**

Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in light, electron, or scanning probe microscopy, or the biological or the physical sciences, MSA takes your knowledge to the next level!

Members Receive:

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and more!



Join MSA Today!

Visit: www.microscopy.org or call 1-800-538-3672



Precision, Speed, Stability

NANO-POSITIONING FOR MICROSCOPY & IMAGING




Manual Stage
25x25mm Travel




High-Speed Nanofocus
<1nm Resolution




XYZ Flexure Stage
<1nm Resolution



UHV Piezo Stage
0.1nm Resolution, XYZ



XY Piezo Motor
Stage, 100mm




The Broadest and Deepest Portfolio of Precision Motion Technologies



Vacuum Stages




Flexure Stages




Air Bearings



Fiber Alignment



Hexapods



Beam Steering

Learn more ► www.pi-usa.us
508.832.3456
PI (Physik Instrumente)

THE NEW DESK V THE GOLD STANDARD IN SAMPLE PREPARATION



Now upgraded with storage for up to ten recipes and a color touch screen display, the Desk V HP delivers breakthrough electron microscopy sample preparation.

DENTON VACUUM



BARRIERS BECOME BREAKTHROUGHS

Visit us at:
www.dentonvacuum.com/mt



SIA

1 TO 50 MEGAPIXELS live and slow scan
MAGNIFICATION FACTOR OF 1 on bottom mounted cameras
DIFFRACTION BEAM STOP on side mounted cameras



Affordable TEM camera systems for research, education, healthcare, and industry since 2001

Scientific Instruments and Applications
2773 Heath Lane • Duluth, GA 30096
(770) 232 7785 • www.sia-cam.com

INDEX TO ADVERTISERS

Bruker Nano GmbH	A5
Denton Vacuum	A15b
Diatome	Cover 2
Electron Microscopy Sciences, EMS	Cover 3
FEI Company	A6
Hitachi	A9
JEOL	Cover 4
Microscopes International, LLC	A13
PI (Physik Instrumente)	A15a
Scientific Instruments & Applications (SIA)	A16a
TESCAN	A11
ZEISS	A3

Introducing the EMS

poly III

the solution for Evaporation-Controlled Automated Embedding and Polymerization

- Automates embedding
- Reduces hands-on time
- Minimizes exposure to hazardous and irritating chemicals
- Reduces solvent and resin use
- Facilitates the processing of up to 52 samples in one instrument run
- Prevents specimen loss

The EMS POLY III is an instrument for the embedding of specimens by the proper combination of pressure and temperature. Central to the instrument is a specimen chamber that is temperature controlled and which can be heated up from room temperature to 70°C. The pressure in the chamber can be reduced from ambient pressure to a controlled level with an inbuilt vacuum-pump. The instrument chamber accepts up to 52 BEEM specimen vials, and features preset programs which can be modified according to the user's preference. In the presets pressure and temperature settings have been coordinated and optimized for an efficient removal of solvent from the specimens. Bulk removal of solvent is followed by steps for the thorough removal of trace amounts.

As a practical approach the instrument can be loaded by the end of a workday and (when using acetone or propylene oxide) by the next morning the vials are ready for polymerization after the vials have been topped up with pure resin.

A lengthy and sometimes tedious manual procedure now reduced to a few simple steps.



Yeast cells were fixed with glutaraldehyde in cacodylate buffer, washed in distilled water and postfixed with 1% KMnO4 in distilled water




EMS has it!

CONTACT US FOR
MORE INFORMATION...


**Electron
Microscopy
Sciences**

P.O. Box 550 • 1560 Industry Rd.
Hatfield, Pa 19440
Tel: (215) 412-8400
Fax: (215) 412-8450
email: sgkcck@aol.com
or stacie@ems-secure.com

www.emsdiasum.com



Facilitates the processing of up to 52 samples in one instrument run.



A choice of 3 embedding programs in the EMS POLY III for 3 different solvents. They have been pre-programmed for general use but the user can change the programs to fit specific specimens.

GAIN FRESH INSIGHT

IMAGING AND ANALYSIS
FROM THE NANO TO
MACRO SCALE

100 μ m

*

100 nm



NEW!
*JSM-IT100 —
up to 300,000X, SE-BSE-EDS,
LV/HV, versatile compact SEM*



NEW!
*JSM-7200F —
10X to 1,000,000X, compact
nanolaboratory, ultrahigh
resolution FE SEM*



NEW!
*JCM-6000Plus NeoScope —
10X-60,000X, SE-BSE-EDS, LV/HV*

**Call us for a demo today!
See our full suite of JEOL SEMs.**

* Courtesy of José R. Almodóvar Rivera, University of Puerto Rico Mayagüez Campus

JEOL
Solutions for Innovation

www.jeolusa.com • salesinfo@jeol.com