

Sub-Micrometer X-ray Tomography of Radiolarians: Computer Modeling and Skeletonization

Roger C. Wagner,^{1*} John R. Jungck,¹ and Denis Van Loo²

¹Department of Biological Sciences, 005 McKinly Lab. University of Delaware, Newark, Delaware 19716

²UGCT-N12, Proeftuinstraat 86, 9000, Ghent, Belgium

*rags@udel.edu

Introduction

X-ray tomography is used to visualize the surface and interior structure of solid objects on the basis of X-ray attenuation (the proportion of X rays scattered or absorbed as they pass through an object) [1]. A tomographic series is acquired by rotating an X-ray source around a specimen and collecting transmitted X-ray intensity on a series of detectors. The detectors acquire a set of images, which are X-ray projections. These are converted to a voxelated 3D file that can be visualized and analyzed with a wide variety of 3D-based imaging programs. Resolution of these image files (designated as voxel pitch) depends on X-ray intensity, focal spot size of the X-ray source, detector size, and noise reduction. Using a combination of these factors, X-Ray Engineering, Belgium XRE (www.xre.be), has developed custom-made X-ray scanning systems designed for a variety of applications. These instruments generate tomographic information well into the sub-micrometer range.

Radiolarians are amoeboid protists with elaborate, radially symmetrical mineral skeletons [2]. Fossilized radiolarian skeletons (also called tests) consist of SiO₂, which is X-ray attenuating. This provides the mechanism for acquisition of 3D data sets by X-ray tomography. With the development of high-resolution computerized tomography (CT), detailed 3D data sets of radiolarians can be acquired with sub-micrometer resolution. We report here the recording of X-ray tomograms (0.103 μm voxel pitch) of selected radiolarian tests and subsequent computer modeling of 3D data sets that reveal intricate details of the exterior and interior of the specimens. Also, digital skeletonization of some of the rendered models reveals their essential geometrical features, which can be measured with great accuracy based on the voxel size of the original data set.

Methods and Materials

Washed radiolarians [3] were obtained from Dr. Richard Howey (University of Wyoming). These specimens originated from sediments collected from the seas of Barbados and were dated to the Tertiary period.

X-Ray tomography. A custom-made, high-resolution CT imaging system, built via a collaboration between UGCT (Ghent University Centre for X-ray Tomography) and its spin-off company, X-Ray Engineering (www.xre.be), was used to acquire tomographic data sets of radiolarians. Several thousand radiolarians were mailed to Ghent, Belgium, on the sticky tab of a post-it-note. Single skeletons or groups of specimens were selected under a stereoscopic microscope and mounted on the tip of a 0.9 mm carbon

sample-holder stick using nitrocellulose as an adhesive. The sticks were inserted into the CT system on a custom sample holder. The X-ray CT system consisted of a transmission X-ray source with a 1 μm thick tungsten target, a high-resolution CCD-based sensor, and eight motorized stages on a damped optical table in a temperature-controlled, shielded room.

The detector was used in a 995×664 pixel mode with a 36 μm pixel size. The X-ray source operated at 50 kV and 0.7 W. For each sample, typically 1,600 two-second projections were acquired over a rotation of 360 degrees. The scan time was approximately 45–60 minutes, using the acquisition software ACQUILA (www.xre.be). The recorded data were processed and reconstructed using XRE dynamic reconstruction code into data sets with voxel size between 300 and 950 nm.

Computer modeling and rendering. 3D image stacks (DICOM or TIFF) were thresholded and surface-rendered with AMIRA 5.6 (FEI, Houston, TX) running on a Mac-Pro platform. Rendered models were rotated and sliced digitally at several angles to reveal the specimen's internal structure. The motion of a virtual camera was scripted to move around and within the modeled radiolarian test to produce a fly-through video for closer assessment of the interior structure.




Figure 1: Rendered model of a group of radiolarians embedded in nitrocellulose adhesive. Several forms are present exhibiting cage-like tests with polygonal elements.

Incomparable...



...and still
innovating

DiATOME diamond knives

ultra 45° • cryo • histo
ultra 35° • histo jumbo
cryo immuno • ultra sonic
ultra AFM & cryo AFM
trimtool 20 • trimtool 45
trimtool 90

Over 40 years of development,
manufacturing, and customer service

DiATOME U.S.

P.O. Box 550 • 1560 Industry Rd. • Hatfield, Pa 19440
Tel: (215) 412-8390 • Fax: (215) 412-8450
email: sgkcck@aol.com • www.emsdiasum.com

Free customer service

Sectioning tests with biological and material research specimens of all kinds. We send you the sections along with the surfaced sample, a report on the results obtained and a recommendation of a suitable knife. Complete discretion when working with proprietary samples.

Re-sharpening and reworking service

A re-sharpened Diatome diamond knife demonstrates the same high quality as a new knife. Even knives purchased in previous years can continue to be re-sharpened. The knives can be reworked into another type of knife for no extra charge, e.g. ultra to cryo or 45° to 35°.

Exchange service

Whenever you exchange a knife we offer you a new DiATOME knife at an advantageous price.



trimtool 90


Many requests from customers doing FIB cutting of biological and technical sample blocks have motivated us to relaunch the trim 90 blade:

Please contact us for more information.

DiATOME manip

Our new tool for an easy handling and mounting of section ribbons!





Figures 2a–2d: Selected radiolarian models exhibiting highly resolved, sharply defined surface elements. Models a, c, and d were volume-edited from the group exhibited in [Figure 1](#).




Figure 3: Digital bisection of rendered models reveal medullary shells within the outer (cortical) shells. In some cases, a single medullary shell was held in place by struts oriented in three planes (3a). In other specimens, two medullary shells were present, each anchored to each other and to the cortical shell by multiple struts (3b).

Surface-rendered models were skeletonized. This is an erosion or thinning of voxels so that only a string of voxels remains. The resulting skeleton represents the essential

geometry of the radiolarian cage-like test. The voxel size of the original data set was used to measure the lengths of the skeleton segments. Segment number and nodes (intersecting points of segments) were also computed. All these functions are available with the latest version of AMIRA (6.0).

The rendered models were converted to VRML files with AMIRA and then to STL files with netfab. Several high-resolution models were 3D printed by Shapeways with a strong, flexible polymer.

Results

Several radiolaria embedded in the nitrocellulose adhesive were scanned as a group ([Figure 1](#)). A variety of forms were present, each perforated with different patterns forming cage-like structures with polygonal outlines. Selected specimens were volume edited (cropped in three dimensions) and surface rendered as individual radiolarians ([Figures 2a–2d](#)). These rendered models exhibited sharply defined features, and their surface details were resolved well into the sub- μm range.

Digital bisection of rendered models revealed medullary shells held in place by struts connected to the outer cortical shell ([Figure 3a](#)). In some specimens, two medullary shells were present, one inside the other. The two inner shells were connected by struts, and both were connected to the outer shell by additional struts ([Figure 3b](#)).

Several of these models were amenable to digital skeletonization [4–6]. This is a thinning algorithm that erodes voxels away from the modeled surface until only a single line of voxels remains (the medial axis). This worked well for specimens with cylindrical or near cylindrical segments, forming the polygons. In such cases the resulting skeleton accurately reflected the essential geometry of the radiolarian test ([Figure 4](#)).

Topological features of the skeletons, such as nodes (intersections of segments) and segment number, were computed. Segment length was calculated according to the voxel size




Figure 4: Digital skeleton (medial axis transform) of a radiolarian test shown against a bisected background of the original model. This represents the essential geometry of the test, including polygonal segments and nodes (intersections of segments). Segment length is derived from the voxel dimension of the data set (inset).

of the file because each segment was a single row of voxels (inset of Figure 4).

Figure 5 shows a 3D-printed model of a radiolarian, large enough to hold in your hand. Currently these models are printed in a strong but flexible polymer.

Discussion

When microscopic biological structures are fluorescent or X-ray attenuating, 3D data sets can be acquired from them. In confocal fluorescence microscopy, a series of confocal fluorescent planes is used to generate a 3D file (z-stack). However, because the exciting fluorescent beam is reduced in energy as it passes through the specimen, the thickness of the z-stack with sufficient fluorescent signal is limited [7]. This limitation is overcome by X-ray tomography, in which




Figure 5: 3D-printed model of a radiolarian produced from a high-resolution tomographic dataset.

a series of X-ray projection images of internal structure are acquired through a rotation of 360 degrees [8].

Micro-CT systems have limited resolution, and sub-micrometer or Nano CT systems are required to resolve fine structure such as corrosion-casted blood capillary systems [8]. Micro-fossils of radiolarians and foraminiferans [9, 10] are ideal for Nano CT because they have X-ray attenuating skeletons with very fine detail, only resolved at the submicron level.

Since the first detailed descriptions of radiolarians by Ernst Haeckel [11], thousands of extant and extinct species have been described [12] based on their collection sites and the geometry of their tests. Early reports relied on microscopy with poor capabilities in contrast, resolution, and depth of field. Scanning electron microscopy has greatly

improved the ability to resolve radiolarian fine structure [13]. Although the SEM can be engineered to view specimens through 360 degrees, the internal structure of radiolarians such as the medullary tests [13] is obscured or hidden from view.

These limitations are overcome with sub-micrometer X-ray tomography, where highly resolved 3D models of radiolarians can be rotated through any angle of viewing and digitally sliced to reveal medullary shells and their associations with outer cortical shells. This greatly improves the precision of radiolarian description and taxonomic classification.

Mathematical analysis based on the geometry of the tests can also be performed because voxel dimension provides a precise measure of polygonal segments of the tests, as well as nodes and segment numbers. The rendered models constitute 3D files, which are amenable to scaling up and used for 3D printing of models that can be held in one's hand.

Conclusion

Sub-micrometer X-ray tomography revealed the intricate structure of siliceous radiolarians at a resolution not achieved by conventional light microscopy. The ability to digitally slice rendered models of radiolarian tests provides additional internal structural details of medullary tests and the struts holding them to the outer wall. Skeletonization of radiolarian models provides a means to quantitate connecting nodes of polygonal units, as well as the number and lengths of the polygonal sides. These capabilities greatly improve structural characterization of radiolarians and their taxonomic classification.

Acknowledgements

We are greatly indebted to Dr. Richard Howey (University of Wyoming) for providing samples of washed radiolarians.

References

- [1] B Flannery et al., *Science* 18 (1987) 1439–44.
- [2] OR Anderson, *Radiolaria*, Springer, New York, 1983.
- [3] A Sanfilippo et al., *Plankton Stratigraphy*, Cambridge University Press, New York, 1985, 107–32.
- [4] S Lobregt et al., *IEEE Transaction on PAMI* 2 (1980) 75–77.
- [5] Y Zhou et al., *The Visual Computer* 14(7) (1998) 303–14.
- [6] K Palagyi et al., *TASK Q* 3 (1999) 397–408.
- [7] R Wagner et al., *Microsc Microanal* 12 (2006) 262–68.
- [8] R Wagner et al., *Microsc Microanal* 17 (2011) 215–19.
- [9] F Mees et al., *Applications of X-ray Computed Tomography in the Geosciences*, Geological Soc Lond Special Pub (2003) 215–30.
- [10] RP Speijer et al., *Geosphere* 21 (2008) 760–63.
- [11] E Haeckel, *Rep Sci Res Voyage H.M.S. Challenger 1873–1876* 18 (1887) 1–1803.
- [12] J Pawlowski and FJ Burke, *Eukaryot Microbiol* 5(1) (2009) 16–25.
- [13] OR Anderson. in *Silicon and Siliceous Structures in Biological Systems*, Springer, New York, 1981, 347–80.

MT




YOUR PARTNER IN X-RAY IMAGING

Custom tailored state-of-the-art X-ray imaging solutions


Modular and flexible multiresolution X-ray CT systems down to nanometer scale resolutions

Dedicated 4D dynamic imaging, bringing synchrotron speeds to the lab

Versatile micro-CT systems optimized for core and long-object scanning




www.XRE.be



bringing correlative light and
electron microscopy forward

CryoCapsule[®]


Designed to accelerate, facilitate and standardize
sample manipulations throughout the CLEM workflow.



EMS has it!

CryoCapsule[®]
Sample Preparation Device
for Correlative Light to
Electron Microscopy

The CryoCapsule[®] is a new tool in the field of High Pressure Freezing (HPF) and correlative light and electron microscopy (CLEM). Comparable to a small petri dish, it is composed of a landmarked sapphire disc and a gold spacer ring (50µm thick) maintained together by a plastic ring^[1].



The specimens are encapsulated between the support sapphire disc (carbon landmarked) and a covering sapphire disc.

The CryoCapsule[®] is loaded into a specific adaptor and live cell imaging can then be done directly on the specimen prior to HPF.

Post-HPF, the specimen is processed for freeze substitution^[2] and room temperature sectioning.

CryoCapCell has also developed a set of tools to manipulate the CryoCapsule[®] [2] in most scientific environments.

**CONTACT US FOR
MORE INFORMATION...**

**Electron
Microscopy
Sciences**

**P.O. Box 550 • 1560 Industry Rd.
Hatfield, Pa 19440**

Tel: (215) 412-8400

Fax: (215) 412-8450

email: sgkcck@aol.com

or stacie@ems-secure.com

www.emsdiasum.com

References

1. Heiligenstein X, Heiligenstein J, Delevoye C, Hurbain I, Bardin S, Paul-Gilloteaux P, Sengmanivong L, Régnier G, Salamero J, Antony C, Raposo G. The CryoCapsule: Simplifying Correlative Light to Electron Microscopy. Traffic [Internet] 2014 [cited 2014 May 14];15:700-16.
2. Heiligenstein X, Hurbain I, Delevoye C, Salamero J, Antony C, Raposo G. Step by step manipulation of the CryoCapsule with HPM high pressure freezers. Methods Cell Biol [Internet] 2014 [cited 2014 Nov 27];124:259-74.

