



Microscopy TODAY

Volume 26 Number 4 2018 July



ETHOS

Focused Ion and Electron Beam System

HITACHI
Inspire the Next

https://doi.org/10.1017/S155192951800055X

...Aspire to Excellence

The All-New Hitachi ETHOS SEM-FIB combines ultra-high resolution imaging and elemental analysis at low voltages with ion optics for nm-scale precision processing.



Designed with a variety of applications in mind


- Large analytical specimen chamber for numerous accessories
- Automated operation including macro processing
- Time sharing mode for dual simultaneous live imaging and processing
- Real-time analytical 3D segmentation capability

Advanced Microsampling

- Sample orientation control with Anti Curtaining Effect (ACE) technologies
- 4-axis lift-out function for advanced TEM specimen preparation

Triple-beam system yielding highest-quality results

- Low acceleration voltage processing with noble gas ion beam
- Selectable ion species (argon/xenon)



200 kv ADF STEM Image of processed lamella by Triple Beam Ar Ion at 1kV

For more information, visit www.hitachi-hightech.com/us/fib



Science for a better tomorrow

Hitachi High Technologies America, Inc.

www.hitachi-hightech.com/us

Tel. 800-253-3053

© 2018 Hitachi High Technologies America, Inc. All rights reserved.

PAPER SUBMISSION PORTAL:

Opens December 2017

SUBMISSION DEADLINE:

February 15, 2018



M&M 2018 **MICROSCOPY &** **MICROANALYSIS**

August 5-9, 2018 • Baltimore, MD



www.microscopy.org/MandM/2018

500 poster presentations over 4 days

>120 companies in the largest microscopy-related exhibit hall in the world

70 FREE Vendor Tutorials

35 named symposia in Physical, Biological, and Analytical Sciences


7 traditional day-long Sunday Short Courses

3 in-depth, pre-meeting Congresses on featured topics, featuring platform and poster presentations and networking

2 fascinating Plenary Session speakers (on 50-cent microscopes, and finding micrometeorites on your roof!)

+ 1 vibrant, historic, waterfront city (with 650 daily flights from 70 destinations)

= M&M 2018 in BALTIMORE!
DON'T MISS OUT!



Aluminum nitride ceramic with yttrium inclusions.

Prisma E SEM for Materials Science

The new Prisma E scanning electron microscope (SEM) combines a wide array of imaging and analytical modalities with new advanced automation to offer the most complete solution of any instrument in its class. It is ideal for industrial R&D, quality control, and failure analysis applications that require high resolution, sample flexibility and an easy-to-use operator interface. Prisma E SEM succeeds the highly successful Quanta SEM.



Find out more at thermofisher.com/EM-Sales

ThermoFisher
SCIENTIFIC

Contents

Light Microscopy

- 12** *Artemia*: A Model Specimen for Educational Microscopy Projects in Biological and Ecological Fields
Jörg Piper
- 20** Guidelines for Understanding Magnification in the Modern Digital Microscope Era
James A. DeRose and Michael Doppler

Electron Microscopy

- 34** Stereoscopic Effects from Single SEM Images
Robert Sturm

Meeting Preview

- 40** IMC19: 19th International Microscopy Congress in Australia
Nicola Thomas

Pioneers of Microscopy

- 44** Cryosectioning and Immunolabeling: The Contributions of Kiyoteru Tokuyasu
Gareth Griffiths, Jan-Willem Slot, and Paul Webster

Highlights from *Microscopy and Microanalysis*

- 56** Assessing Soft-Tissue Shrinkage Estimates in Museum Specimens Imaged With Diffusible Iodine-based Contrast-Enhanced Computed Tomography (diceCT)
- 56** Characterization of Lithium Ion Battery Materials with Valence Electron Energy Loss Spectroscopy
- 57** Electron source brightness and illumination semi-angle distribution measurement in a transmission electron microscope

About the Cover



Morphological differentiation of *Artemia franciscana* during development from birth to full adult. Top: day 1 (left) and day 13 (right), bottom: day 22 (left) and day 35 (right). Various magnifications.

See article by Piper.

Departments

- | | |
|--------------------------------------|--------------------------------|
| 7 Editorial | 60 NetNotes |
| 8 Carmichael's Concise Review | 68 Calendar of Meetings |
| 50 Industry News | 72 Dear Abbe |
| 52 Product News | 74 Index of Advertisers |



Providing microscopy supplies and specimen preparation equipment to our valued customers for half a century.

Visit us August 3–5 in Baltimore, MD
M&M 2018 | BOOTH #1012



NEW SBEM Tools & Supplies
FEI VolumeScope Pin Stub, Storage, Tweezers & Colloidal Silver Paint




Cryo-EM Tools & Accessories
Large Dewars, Foam Dewars, Grid Boxes & Grippers



PELCO® Modular SEM Holders & Mounts



NEW Products for Specimen Thinning in Preparation for Electron Microscopy




NEW Quantifoil® TEM Substrates



Cressington Coating Systems
Carbon Evaporation & Metal Sputtering




PELCO BioWave® Pro+ Microwave Tissue Processor




M&M 2018 Vendor Tutorial

Join us Tuesday, August 7 at 5:45pm, Booth #1012



PELCO BioWave® Pro+ Assisted Serial Block-Face SEM (SBEM)

Using the PELCO BioWave® Pro+, the entire SBEM process can now be completed in under a day with identical results to the standard long protocol.



Images: Rick Webb, University of Queensland



WIN \$50 OFF YOUR NEXT ORDER!

Visit us at M&M 2018, and enter to win one of fifty gift certificates for \$50 off your next order.



www.tedpella.com sales@tedpella.com 800-237-3526