





# Instant Results for Your Life Science TEM Lab

GATAN  
TEM  
DIGITAL  
IMAGING



Tomorrow's  
Technology  
Today



## Instant Results from Gatan CCD Cameras

Gatan offers a wide range of CCD cameras and the most productivity software for the TEM microscopist. Based on our "industry standard" DigitalMicrograph™ image acquisition, manipulation, and archiving platform, the microscopists have many expansion options, including

- Automatic TEM alignment, stigmation, and focus adjustment
- Digital Montage display
- Automated LowDose operations

*Gatan offers a wide range of hardware and software for electron microscopy, including a wide range of CCD cameras, electron energy loss imaging and spectroscopy systems, specimen holders and prep equipment. Visit our web site or contact a local representative for more information.*





**Gatan Inc. (USA)**  
Tel 925 463 0200  
Fax 925 463 0204  
**Gatan Inc. (USA)**  
Tel 724 776 5260  
Fax 724 776 3360

**Gatan Online**  
<http://www.gatan.com>  
[info@gatan.com](mailto:info@gatan.com)  
[help@gatan.com](mailto:help@gatan.com)

**Gatan GmbH (Germany)**  
Tel 089 352 374  
Fax 089 359 1642  
**Gatan Ltd. (UK)**  
Tel 01536 743150  
Fax 01536 743154

**Gatan Inc. (Singapore)**  
Tel: 65 293 3160  
Fax: 65 293 3307  
**Nippon Gatan (Japan)**  
Tel: 81 03 5387 3585  
Fax: 81 03 5387 3585

# As Easy As It Looks



PEO at 60°C with evidence of melting clearly shown after 2 minutes.

- ▶ **Based on a Simple PC Interface.**
- ▶ **Intuitive Operation with a Mouse Alone.**
- ▶ **Capable of Imaging Under Either Ambient or High Vacuum and at Both High and Low Temperatures.**
- ▶ **Fully Vibrationally Isolated.**
- ▶ **Capable of Wet or Dry Imaging.**
- ▶ **Variable Gas Pressure/Temperature.**

The new JSPM-4200 is a multipurpose turnkey Scanning Probe Microscope (SPM) capable of imaging using a myriad of different modes and

JEOL USA, Inc., 11 Dearborn Road, Peabody, MA 01960  
Tel: 978/535-5900 • Fax: 978/536-2205  
Email: [spm@jeol.com](mailto:spm@jeol.com) • WEB: <http://www.jeol.com>

## NEW - Multi-Environmental

# SPM

## Scanning Probe Microscope

accessories. Imaging modes include STM, STS (Scanning Tunneling Spectroscopy), CITS (Current Imaging Tunneling Microscopy), AFM, LFM (Lateral Force Microscopy), Non-contact Mode AFM, Discrete Contact Mode AFM, Lateral Modulation LFM, MFM (Magnetic Force Microscopy), and Force Modulation Microscopy.

If you would like to learn more about JEOL SPM's please give us a call and arrange for your own hands-on demonstration so that we can prove to you that it really is "As Easy As It looks".

# JEOL