

COMMITTED TO EXCELLENCE

Perfectly smooth sections
from 25nm to 5 μ m.

Hydrophilic diamond
surface draws water to the
edge for easy sectioning.

Flawless, jewel quality
diamond crystal, in sizes
from 1 to 10mm.

Maximum sectioning stability is
assured by molecular welding
between diamond and stainless
steel shank.

Highly inert and strong bonding
material is unharmed by solvents.



Nine different boat styles and seven types
for all ultramicrotomy, cryo, histology
and materials applications.

MICRO STAR DIAMOND KNIVES

Tel: 800 533 2509 Fax: 409 294 9861

E-mail: mistar@msn.com

For full prices and information visit our web site: <http://www.microstartech.com/>

Circle Reader Inquiry #25