

The first of a series of announcements describing electro-surgical technique with Marconi Diathermy equipments as the operative medium.

Transverse Electro-Surgical Incision

- 1** Cutting needle — preferably tungsten or fine steel.
- 2** Note angle of cutting needle.
- 3** Note absence of excessive sparking at point of contact.
- 4** Wound edges dry and covered with fine film of coagulum.
- 5** First incisional depth does not extend below subcutaneous fat.

Marconi Instruments Ltd. now have in production a wide range of Diathermy electro-surgical equipments designed for Ophthalmic, Brain, Urological and General Surgery. Full particulars, together with information upon their applications, are freely available on request.


MARCONI  **INSTRUMENTS LTD.**

Electro-Medical Apparatus

ST. ALBANS, HERTFORDSHIRE · TELEPHONE: ST. ALBANS 4323/6

Northern Office: 30 ALBION STREET, HULL · TELEPHONE: HULL 16144
 Western Office: 10 PORTVIEW ROAD, AVONMOUTH, BRISTOL · TELEPHONE: AVONMOUTH 43R
 Southern Office: 109 EATON SQUARE, LONDON, S.W.1 · TELEPHONE: SLOANE 8615

Please mention *The Journal of Laryngology* when replying to advertisements



Instruments used in the practice of Laryngology, Rhinology
and Otology

Mayer & Phelps

MANUFACTURERS
of Surgical Instruments
and Aseptic Hospital
Furniture

Established over
Eighty Years

CHIRON HOUSE

32-34 New Cavendish St., London, W.1

DOWN BROS. and MAYER & PHELPS have amalgamated
The personal managements remain as heretofore

DOWN BROS. and MAYER & PHELPS LTD.