Compositional and Chemical Segregation in $Li_{1.2}Ni_{0.2}Mn_{0.6}O_2$ Cathode Materials Characterized by Atom Probe Tomography and Scanning Transmission X-ray Microscopy

A. Devaraj¹, R. Colby¹, M. Gu¹, C.M. Wang¹, T. Tyliszczak⁴, J. Zheng², J. Xiao², J. Zhang², I. Belharouak³, D. Wang³, K. Amine³, S. Thevuthasan¹

- 1 Environmental and Molecular Sciences Laboratory, Pacific Northwest National Laboratory, Richland WA, 99354
- 2 Energy and Environmental Directorate, Pacific Northwest National Laboratory, Richland WA 99354
- 3 Chemical Sciences and Engineering Division, Argonne National Laboratory, Argonne IL, 60439
- 4 Advanced Light Source, Lawrence Berkeley National Laboratory, Berkeley CA, 94720

Li-ion batteries are one of the most widely used rechargeable energy sources, with a two times higher voltage than aqueous batteries, six times higher gravimetric specific energy compared to the lead acid battery [1], and a high energy density [2]. Layered Li_{1.2}Ni_{0.2}Mn_{0.6}O₂ (LMNO) nanoparticles have shown to be a very promising next generation high capacity cathode material for Li-ion batteries. However, there is still a debate as to the structure, cation ordering and distribution in the cathode, both before and after cycling [3,4]. In addition, the layered cathode has shown a trend to transform to spinel phase after high voltage cycling [5]. This important structural change would result in capacity fading and poor rate performance, limiting its application in heavy-duty electric vehicles. These important structural, compositional and chemical changes can be clarified using a combination of atom probe tomography (APT), scanning transmission x-ray microscopy (SXTM) and transmission electron microscopy (TEM) analysis.

Recent scanning TEM x-ray energy dispersive spectroscopy (STEM-EDS) tomography results have indicated the segregation of Ni at surfaces and grain boundaries. It has been proposed that this clustering acts as a barrier to the Li [6]. Moving forwards, an effective means of detecting Li is needed. Laser-assisted, local electrode atom probe (LEAP) combines a time of flight mass spectrometer with a point projection microscope capable of position sensing to perform atomic scale imaging of materials with sub-nanometer (<0.2 nm) spatial resolution in three dimensions and parts-per-million (ppm) composition sensitivity [10], and is capable of detecting Li. LMNO nanoparticles have been prepared for APT analysis by a focused ion beam liftout and analyzed in a LEAP 4000XHR with a UV laser. In addition to confirming Ni and Mn segregation observed previously be STEM-EDS, the atom probe reconstructions demonstrate that the Li preferentially segregates to Mn-rich regions over the Ni-rich regions (Fig. 1a,b). In addition, LMNO before and after cycling were examined using STXM on Advanced Light Source beamlines 5.3.2.1 and 11.0.2. Examination of the oxygen K-edge evidenced the emergence of a peak at 533.5 eV in samples after cycling, predominantly located at surfaces. This feature is additionally observed in spinel LiMn_{1.5}Ni_{0.5}O₄ materials only after cycling, and has been previously linked to lithiation, potentially suggesting a surface-rich segregation [7,8].

- [1] P. Verma, P. Maire, and P. Novak, *Electrochimica Acta*, **55** (22), 6332 (2010).
- [2] J. Liu and A. Manthiram, Journal of Materials Chemistry, 20 (19), 3961 (2010).

- [3] Jarvis, K.A., et al.. Chem. Mater., 2011. 23, 3614-3621.
- [4] Armstrong, A.R., et al., Journal of the American Chemical Society, 2006. 128, 8694-8698.
- [5] Hong, J., et al., J. Mater. Chem., 2010. 20, 10179-10186.
- [6] M. Gu, I. Belharouak, A. Genc, Z. Wang, D. Wang, K. Amine, F. Gao, G. Zhou, S. Thevuthasan, D. R. Baer, J. Zhang, N.D. Browning, J. Liu, and C.M. Wang. *NanoLetters* Sept. 17, 2012
- [7] A portion of the research was performed using EMSL, a national scientific user facility sponsored by the DOE-BER and located at PNNL.
- [8] The Advanced Light Source is supported by the Director, Office of Science, Office of Basic Energy Sciences, of the U.S. Department of Energy under Contract No. DE-AC02-05CH11231.




Figure 1: (a) APT reconstructions separately showing the Li, Mn, and Ni ions for a typical unicycled Li_{1.2}Ni_{0.2}Mn_{0.6}O₂ particle and (b) density maps from a sliced subsection of the reconstruction demonstrating the segregation of the Mn and Ni, along with correlation of the relative density of the Li with the Mn. (c) A three-color map extracted from a STXM stack of the oxygen K edge for cycled LMNO, demonstrating the spatial variation of features at ~529.2 eV (red), ~531.6 eV (blue), and 533.5 eV (green). The scale bar is 600 nm.