




**VACUTAINER®
QuikHeel™
Safety Lancet**



**VACUTAINER®
Safety-Lok™
Needle Holder**



**BD® Saf-T E-Z Set®
Winged Infusion Set**



**BD® Safety-Lok™
Syringes**



**BP® Bard-Parker®
SAFETYLOCK™
Surgical Blade System**



**VACUTAINER®
SafetyGlide™ Blood
Collection Assembly**




**BD® Saf-T E-Z Set®
Winged Infusion Set**




**BD® Cannula
InterLink®
IV Access System**



**VACUTAINER®
Safety-Lok™
Needle Holder**




**BD® Insite®
AutoGuard™
Shielded IV Catheter**



**BD® Safety-Lok™
Syringes**



**BD® GUARDIAN®
Sharps Collector**



© 1999 Becton Dickinson and Co. For more information, call 1-888-237-2762, press Prompt #1 or visit our web-site at <http://www.bd.com>

For more information, call 1-888-237-2762, press Prompt #1 or visit our web-site at <http://www.bd.com>


Insite is a trademark of Baxter International Inc.
All other noted trademarks are property of Becton Dickinson and Co., Franklin Lakes, NJ

© 1999 Becton Dickinson and Co.

Your safety. Our concern.

Becton Dickinson, an innovator of quality medical products for more than a century, is also the leader in sharps safety. Our products, educational programs and exposure surveillance systems help to provide you with a safer work environment, whether you are giving injections, drawing blood, or disposing of sharps.

BD. One Company, Many Solutions.




**BD® Insite®
AutoGuard™
Shielded IV Catheter**




**BD® SafetyGlide™
Shielding
Injection Needle**




**BD® GUARDIAN®
Sharps Collector**



**BD® Saf-T-Intima™
IV Gatherer
Safety System**



**BD® Cannula
InterLink®
IV Access System**



**VACUTAINER®
Safety-Lok™
Blood Collection Set**



For more information, call 1-888-237-2762, press Prompt #1 or visit our web-site at <http://www.bd.com>

InterLink is a trademark of Baxter International Inc.
All other noted trademarks are property of Becton Dickinson and Co., Franklin Lakes, NJ
<https://doi.org/10.1017/S0195941700068764> Published online by Cambridge University Press

© 1999 Becton Dickinson and Co.

NEW in endoscope reprocessing

Greater FLEXIBILITY

We've added a new twist to our hydrogen peroxide solution



SPOROX II Solution—now with a wider range of benefits for greater flexibility

- Improved compatibility with heat-sensitive medical devices containing brass and copper¹
- Significantly less foaming²

Plus all the benefits you love

- Fast-acting sporicide, tuberculocide, fungicide, virucide, and bactericide¹
 - High-level disinfection in 30 minutes at room temperature
 - Sterilization in 6 hours at room temperature
- Does not contain glutaraldehyde
 - No aldehyde odor
- Unique, ready-to-use hydrogen peroxide formula for automatic and manual reprocessing¹
- Oxidative formula reacts with organic matter and facilitates its removal

References:

1. SPOROX II Solution package insert.
2. Data on file, Reckitt & Colman Inc.

SPOROX II
STERILIZING & DISINFECTING SOLUTION
From the makers of LYXOL® I.C.™ products

The simple solution for reprocessing



RECKITT & COLMAN
Professional

The name to trust for your infection control needs

© 1999 Reckitt & Colman Inc.